

Product datasheet for **MC207623**

Jam2 (NM_023844) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Jam2 (NM_023844) Mouse Untagged Clone
Tag: Tag Free
Symbol: Jam2
Synonyms: 1110002N23Rik; 2410030G21Rik; 2410167M24Rik; JAM-2; JAM-B; Jcam2; VE-JAM
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207623 representing NM_023844
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGAGGAGCCCCAAGGCCTCTGATGCTGCTGCTGCTACTACTTGGATCGTCGCCCTGGACTATC
ATAAGGCCAAATGGGTTTTCTGCATCAAAGACCACCGTCAAGAAGTCACAGTAATAGAGTTCCAAGAGGC
TATTTTGGCTTGTAAAACCCCAAAGAAGACTACCTCCTCCAGACTGGAGTGAAGAAGGTGGACAGGGG
GTCTCCTTGGTCTACTACCAACAGGCTCTCAAGGTGACTTTAAAGACCGTGTGAGATGATAGATTCA
ATATACGAATCAAAAATGTTACAAGAAGTGATGCTGGAGAGTATCGCTGTGAAGTCAGCGCTCCGACTGA
GCAAGGCCAGAACCTGCAGGAAGATAAAGTCATGCTAGAAGTACTAGTGGCTCCTGCTGTTCTGCTGT
GAAGTGCCCACTTCTGTTATGACTGGAAGTGTGGTGGAGCTACGATGCCAGGATAAAGAAGGAAACCCAG
CTCCGGAGTACATCTGGTTAAAGATGGCACAAGTTTGTAGGGAATCCAAAAGGCGGCACACACAACAA
CAGCTCGTACACAATGAACACGAAGTCTGGAATTCTGCAATTCAACATGATTTCCAAGATGGACAGTGGA
GAGTATTACTGCGAAGCCCGGAACCTGTGCGGACACCGCAGGTGCCCTGGGAAGCGAATGCAAGTAGATG
TTCTCAACATAAGCGGCATCATAGCAACGGTTGTGGTGGTGGCCCTCGTGATTTCTGTATGTGGCCTTGG
CACATGCTATGCTCAGAGGAAAGGCTACTTTTCAAAGAAACTTCTTCCAGAAGGGCAGTCCTGCATCT
AAAGTCACTACGATGAGCGAAAATGATTTCAAGCACACAAAATCCTTTATAATTTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_023844
Insert Size: 897 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_023844.5](#), [NP_076333.3](#)

RefSeq Size: 4822 bp

RefSeq ORF: 897 bp

Locus ID: 67374

UniProt ID: [Q9J159](#)

Cytogenetics: 16 C3.3

Gene Summary: May play a role in the processes of lymphocyte homing to secondary lymphoid organs. [UniProtKB/Swiss-Prot Function]