

## Product datasheet for **MC207627**

### **Ikkip (NM\_026166) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ikkip (NM\_026166) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ikkip  
**Synonyms:** 120009F10Rik; 1700023M03Rik; D18354; Ikip  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207627 representing NM\_026166  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCCGAAGTGAAGAGCCGAAAGAAGCCGGGGCCAAAGTGGCTGCTCCTGAGCCCGAGAAGCGGAGTG  
ATGGCAGGAAAAACCCGAGGCGCGGGGGACCGGGCTGGGCGGACCCCGCACTGGGCTCAGCCTGCT  
GTCTCTGGCGATGACCCTGGGCCTGGCCTGGCTTGTATTTCAACAATCAGAAAAGTTTGGCAAGGTAGAA  
AAACAATACCGGTTACTGCAAACAGAGAGCAGCGAGTTTCAGGGCCTCCAGAGTAAAATCAGTTTAAATTT  
CCAGCAAGCTTGAGTCTACTGAAAATACCTGCAGGAAGCTACATCATCCATATCTTTGATGACCCAATT  
TGAACAGGAAGTATCTGGCCTCCAAGATCCATACGTGATATTGAGACTAGCGAAGAGATGCTCACACAA  
AAGATGCAAAACCTTAATGAGAAATCCAGAACATTACAGACTTCTGGAAGAGAACAACACTGGCAGAAATGA  
TCGATGACACAGCCGCTTCAAATCAGAAGTGAAGGATACACATTCTGAAGTTACCCTGAAGATTAACCTC  
CGCGGATCAAGAAATAAAATCGCTCACTGAACGTCTAAAAGATTTAGAAGACAGCACACTGCGAAACATC  
AGAACAGTGAGCAGGCAAGAAGAGGAAGATCTTCTACGAGTCGAGGCCAGCTTAGCTGTATACCAAAG  
CAGTTAAGAAGCTAGAAGAAGAACAGCACACGCTCCTAGCCAGAGACGAAGATCTGACCAACAAACTCTC  
CAGCTACGAACCCAAAGTAGAAGAGTGCAAGGCACACTTGCCAACCATAGAAAATGCTGTCCACTCGGTT  
CTCAGAGTCTCTCAGGATCTGATAGGGACAGAAAGGAAGATGGAGGAGCTGACGATGCAGATGTTCAACA  
TGGAGGACGACATGCTGAGGGCTGTGTCTGAGATCATGGAGATGCAGAAAACCTGGAGGGGATCCAGTA  
TGACAACAGCTTACTAAAGATGCAGAACGAGCTGGTCGTTCTAAAAGGAAAAGTTCATGACTTCATAGCA  
TATTCAGTGACCGGAAAAGGGGACTTTGGGAGAGTATAGCCTAGGAAATAAGGGAACCTGATGAGTATT  
AA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_026166
<b>Insert Size:</b>	1122 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_026166.2</a></u> , <u><a href="#">NP_080442.2</a></u>
<b>RefSeq Size:</b>	3251 bp
<b>RefSeq ORF:</b>	1122 bp
<b>Locus ID:</b>	67454
<b>Cytogenetics:</b>	10 C2
<b>Gene Summary:</b>	Target of p53/TP53 with pro-apoptotic function.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).