

Product datasheet for **MC207631**

Hikeshi (NM_026304) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hikeshi (NM_026304) Mouse Untagged Clone
Tag: Tag Free
Symbol: Hikeshi
Synonyms: 0610007P06Rik; 1110002N09Rik; l(7)6Rn; l7Rn6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207631 representing NM_026304
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTTGGCTGCCTGGTGGCGGGAGGCTGGTGCAAACAGCTGCACAGCAAGTGGCAGAGGATAAATTTG
TTTTTGACTTGCCTGATTATGAAAATATCAACCATGTTGTGGTTTTATGCTGGGAACAATCCCATTTC
TGAGGGCATGGGAGGATCTGTCTACTTTTCCTATCCTGATTCAAATGGGGTGCCAGTGTGGCAGCTCCTA
GGATTTGTACGAATGGAAAGCCAAGTGCCATCTTCAAATATCAGGTCTTAAATCTGGGAAGGAAGCC
AGCACCCATTTGGAGCCATGAATATTGTGCGAACCCCATCTGTTGCCAGATTGGCATCTCGGTGGAATT
GTTGGACAGTCTGGCTCAGCAGACTCCTGTAGGCAGTGCTGCCGTGTCCTCCGTTGACTCCTTACGCAG
TTCACACAGAAGATGTTGGACAACCTTACAATTTGCTTCATCATTTGCTCTCTCAGGCCAGATGA
CACCAATCCATCAGAAATGTTTCATCCCAGCAAATGTGGTTCTGAAATGGATGAAAACCTTCAAAGACG
ACTAGCACAGAACCCCTCTTTTGGAAAACA**TAA**

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026304
Insert Size: 594 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026304.3, NP_080580.1</u>
RefSeq Size:	1085 bp
RefSeq ORF:	594 bp
Locus ID:	67669
UniProt ID:	<u>Q9DD02</u>
Cytogenetics:	7 50.4 cM
Gene Summary:	Acts as a specific nuclear import carrier for HSP70 proteins following heat-shock stress: acts by mediating the nucleoporin-dependent translocation of ATP-bound HSP70 proteins into the nucleus. HSP70 proteins import is required to protect cells from heat shock damages. Does not translocate ADP-bound HSP70 proteins into the nucleus (By similarity). May also be indirectly required for organization and/or function of the secretory apparatus in Clara cells in lung.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (a).