

Product datasheet for **MC207696**

Ttc9c (NM_027412) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ttc9c (NM_027412) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ttc9c
Synonyms: 2210019E14Rik; 6330408J23Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207696 representing NM_027412
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGAAGCGCCTTCAGGAGGCTCAGGTATACAAGGAGGAAGGAAACCAGCGCTACCGGAAGGAAAGT
ACCGAGATGCTGTGAGTAGGTACCATCGAGCTCTTCTTCAGCTGCGGGGTCTGGACCAAGTCTGCCTTC
CCCTCTTCTAGTCTGGACCTCAGGGCCAGCTCTCACACCGGAGCAAGAAAACATACTGCACACCATC
CAGACCCACTGTACAACAACCTAGCTGCCTGTCTCCTCCAGATGGAGCCAGTGAACATGAACGAGTGA
GGGAATACAGTCAGAAAGTCCTAGAACGACAACCCGATAATGCCAAGGCCTTGTATCGGGCTGGAGTAGC
CTTCTCCACCTGCAAGACTATGACCGGGCTCGACACCACCTGCTTGCTGCTGCTCAATAGGCAACCTAAA
GATGCCAATGTTTCGGCGGTACCTCCAGTTGACGCAGTCAGAAGTACAGCTACCATAGGAAGGAGAAAC
AGCTGTACCTGGGCATGTTTGGT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_027412
Insert Size: 516 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_027412.3</u> , <u>NP_081688.2</u>
RefSeq Size:	3470 bp
RefSeq ORF:	516 bp
Locus ID:	70387
UniProt ID:	<u>Q810A3</u>
Cytogenetics:	19 A