

## Product datasheet for **MC207698**

### Gfod2 (NM\_027469) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gfod2 (NM\_027469) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gfod2  
**Synonyms:** 5730466C23Rik  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207698 representing NM\_027469  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAAGCTGCTCCAGGAGTCGGTGTGTTCCGGACTGGCAGCTCAGCCCGGTCTGGTCCCGCTGCTGA  
GGGCAGAGGGCTTACCCTGGAGGCCCTGTGGGGAAGACTGAGGAAGAGCAAAGCAGCTGGCCGAGGA  
GATGAACATCACCTTCTACACCAGCCGCACAGACGACGCTTTGCTGCATCAAGACGTAGACTGGTCTGC  
ATCAACATCCACCTCCGCTCACCCGGCAGATATCCGTGAAGGCTCTAGGAATTGGGAAGAACGTGGTGT  
GTGAGAAGGCAGCGACATCCATGGACGCCTTCCGGATGGTGACCCGATCTCGCTACTACCCACAGCTGAT  
GAGCTTGGTGGGGAACGTTCTGCGCTTCTGCCTGCCTTCGTGCGCATGAAACAGCTCATTGCAGAGCAC  
TACGTGGGTGCCGTGATGATCTGCGACGCTCGCATCTACTCCGGGAGCCTGCTCAGCCCCAGCTACGGCT  
GGATCTGTGACGAGCTTATGGGGGGCGGGGCCCTGCACACCATGGGTACCTATATCGTGGACCTGTTGAC  
CCACTTGACTGGCCAGAAAGCCGAGAAGGTGCACGGGCTCCTCAAGACCTTTGTGAGGCAGAATGCCACC  
ATCCGTGGCATCCGCCATGTCACCAGTGATGACTTCTGTTTTTCCAGATGCTCATGGGTGGGGGGTAT  
GCAGCACAGTGACCTCAACTTCAACATGCCAGTGCCTTTGTGCATGAGGTCATGGTAGTGGGCTCCGC  
AGGGCGCCTTGTGCCCGGGGAGCTGACCTCTATGGGCAGAAGAACTCTGCTGCACAGGAGGAGCTGCTG  
GTGAGAGACTCTCTGGTGTGGGCGCAGGGCTCCCTGAGCAAGGTCCCCAGGATGTCCCACTGCTTTACT  
TGAAAGGTATGGTCTATATGGTTCAGGCCTTGCCTCAGTCCTTCCAGGGCAGGGGGATCGCCGTACCTG  
GGACCGCACCCCTGTCTCTATGGCTGCCTCATTGAGGACGGACTGTATATGCAGAGCGTGGTAGATGCT  
ATTTAAAGGTCCAGCCGGTCTGGGGAGTGGGAGACTGTGGAGATGCTGGCAGAGGAGCCAGATGCCAACC  
AGAATCTAAGTGAACACTTCAGCGGAACAACCTG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_027469
<b>Insert Size:</b>	1158 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_027469.4</a></u> , <u><a href="#">NP_081745.1</a></u>
<b>RefSeq Size:</b>	2036 bp
<b>RefSeq ORF:</b>	1158 bp
<b>Locus ID:</b>	70575
<b>UniProt ID:</b>	<u><a href="#">Q9CYH5</a></u>
<b>Cytogenetics:</b>	8 D3
<b>Gene Summary:</b>	Promotes matrix assembly.[UniProtKB/Swiss-Prot Function]