

## Product datasheet for MC207702

### Gpr173 (NM\_027543) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gpr173 (NM\_027543) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gpr173  
**Synonyms:** 3230401K02Rik; Sreb3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207702 representing NM\_027543  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGCCAAACACCACCGGAGAGCCCGAAGAGGTGAGCGGCGCACTGTCTCTGCCATCAGCATCGGCTTATG  
 TGAAGCTGGTGTCTGGGACTGATCATGTGTGAAGCCTGGCAGGCAATGCCATCTTGCCCTGTGGT  
 GCTCAAGGAGCGTGCCCTGCACAAGGCTCCTTACTACTTTCTGTGGACCTGTGCCTAGCCGATGGCATA  
 CGCTCTGCCATCTGCTTCCCCTTTGACTGGCTTCTGTGCGCCATGGCTCCTCATGGACCTTCAGTGCAC  
 TCAGCTGTAAGATTGTGGCCTTTATGGCTGTGCTCTTTGCTTCCATGCGGCCTTCATGCTGTTCTGCAT  
 CAGCGTCAACCGCTACATGGCCATCGCCACCACCGCTTCTATGCCAAGCGCATGACACTCTGGACATGC  
 GCAGCTGTCTGATGGCCTGGACCTTGTCTGTGGCCATGGCTTCCCACCTGTCTTTGATGTGGGCA  
 CCTACAAGTTTATCCGAGAGGAGGACAGTGCATCTTTGAGCATCGCTACTTCAAAGCAAATGACACACT  
 GGGCTTTATGCTTATGTTGGCTGTGCTCATGGCAGCCACACATGCTGTCTATGGCAAGCTGCTACTCTTC  
 GAGTATCGTCACCGCAAGATGAAGCCAGTGCAGATGGTGCCAGCCATCAGCCAAAATGGACATTCATG  
 GCCCTGGGGTACCGGCCAGGCTGCTGCCAAGTGGATCGCTGGTTTTGGCCGTGGGCCCATGCCACCAAC  
 TCTGCTGGGTATCCGGCAGAATGGGCATGCAGCTAGCCGGCGGCTACTGGGCATGGACGAGGTCAAGGGT  
 GAAAAGCAGCTGGGCCGAATGTTCTACGCGATTACACTGCTCTTCTGCTCCTCTGGTCAACCATACATTG  
 TGGCCTGCTACTGGCGAGTGTGTTGTAAAGCCTGCGCTGTGCCACCCTACTTGGCCACTGCTGTTTG  
 GATGAGCTTCCGCCAGGCTGCCGTCAACCCAATCGTCTGCTTCTGCTTAAACAAGGACCTCAAGAAGTGC  
 CTGAGGACTCATGCCCTTGCTGGGGCACAGGAGGTGCCACAGTCCCAGAGAACCTACTGTGTCATGT  
 GA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_027543
<b>Insert Size:</b>	1122 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_027543.4</a> , <a href="#">NP_081819.2</a>
<b>RefSeq Size:</b>	4868 bp
<b>RefSeq ORF:</b>	1122 bp
<b>Locus ID:</b>	70771
<b>UniProt ID:</b>	<a href="#">Q6PI62</a>
<b>Cytogenetics:</b>	X F3
<b>Gene Summary:</b>	Orphan receptor.[UniProtKB/Swiss-Prot Function]