

## Product datasheet for MC207703

### Prss22 (NM\_133731) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Prss22 (NM\_133731) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Prss22  
**Synonyms:** 4733401N09Rik; BSSP-4; SP001LA  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207703 representing NM\_133731  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGATGATCTCCAGACCTCCCCAGCACTGGGTGGGACCAGTTCAGCATCTTAATCCTTCTGGTGCTGC  
 TGACTTCCACAGCTCCCATCAGTGCTGCCACCATCCGAGTGTCCCAGACTGTGGGAAGCCTCAGCAGCT  
 GAACCGGATTGTGGGAGGTGAGGACAGCATGGATGCCAGTGGCCCTGGATTGTTAGCATCTCAAGAAT  
 GGCTCCCACCACTGTGCAGGCTCCCTGCTCACCACCGCTGGGTGGTACAGCCGCGCACTGCTTTAAGA  
 GCAATATGGACAAACCATCTCTGTTCTCAGTATTGTTGGGGCCCTGGAAGCTGGGAGCCAGGCCAAAG  
 GTCCCAGAAAGTAGGCATTGCTTGGGTGCTGCCTCACCCAGGTATTCTTGAAGGAGGGAACCCATGCA  
 GACATTGCCCTGGTGCCTGGAACACTCCATCCAGTTCTCTGAGCGGATCCTGCCATCTGCCTACCTG  
 ACTCCTCTGTCCGTCTCCCTCCAAGACCGACTGCTGGATTGCCGGCTGGGAAGCATCCAGGATGGAGT  
 GCCCTGCCACCCTCAGACCCTCAGAAGCTGAAGGTGCCATCATCGACTCCGAACCTGCAAAAGC  
 TTGTACTGGCGGGAGCCGGTCAGGAAGCCATCACCGAGGGCATGCTGTGTGCTGGTTACCTGGAAGGGG  
 AGCGGGATGCTTGTCTGGGCGACTCTGGGGTCCCTGATGTGCCAGTGGATGACCACCTGGCTACTGAC  
 TGGCATAATCAGCTGGGAGAGGGCTGCGCGGAGCGCAACCGCCGGTGTGTACACCAGCCTCTAGCT  
 CACCGCTCCTGGGTGCAAGGATCGTTCAAGGGTGCAGCTGCGCGGTACTTGGCGGACAGTGGGACA  
 CAGGAAGCTCCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_133731  
**Insert Size:** 924 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_133731.2</a> , <a href="#">NP_598492.2</a>
<b>RefSeq Size:</b>	1343 bp
<b>RefSeq ORF:</b>	924 bp
<b>Locus ID:</b>	70835
<b>UniProt ID:</b>	<a href="#">Q9ER10</a>
<b>Cytogenetics:</b>	17 12.18 cM