

## Product datasheet for **MC207725**

### 2810006K23Rik (NM\_001134717) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 2810006K23Rik (NM\_001134717) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** 2810006K23Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC207725 representing NM\_001134717  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTCTTCTAGGAGCACCTGGGCTCTGCTCCGCCTCCCTCTACCACTGATCAGGATATGCTCTGGGAAGT  
 GGGGGCTCAGGCTTCAGGAGAAGCCAGCACTGCTTCCCAGGAATGGCTGCCAGCACAGTACAGGTGGC  
 AGGCAGGAAGGACTACCTGTCTGTCTCCCTGAATGAGAGTGAGCTCGAAGAACAGTTCGTGAAAGGA  
 CATGGCCAGGGGCCAGGCCACCAACAAGACCAGCAATTGTGTAGTGCTCAAACACGTGCCCTCCGGCA  
 TTGTGGTCAAGTGCCACCAACAAGATCTGTGGATCAAACAGGAAGATAGCTCGGAAAGTCCCTCAGGA  
 GAAAGTGATGTTTTCTACAATGGTGAAAACAGCCCGTTCACAAAGAGAAGCTCGAGGCTGAGAGGAGA  
 AAGCGAGAGAGGAAGAAAAGAGCAAAGGAGACTCTAGAAAAGAAAAAATTATTGAAAGAATTACGGGAAG  
 CAAGTCAAACATCACAGAGAAAAGGCTGATGCTGATGGTATACCAAGGGGCTTCCAGGA**AG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001134717  
**Insert Size:** 555 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001134717.2, NP_001128189.1</u>
<b>RefSeq Size:</b>	1694 bp
<b>RefSeq ORF:</b>	555 bp
<b>Locus ID:</b>	72650
<b>UniProt ID:</b>	<u>Q80VP5</u>
<b>Cytogenetics:</b>	5 F
<b>Gene Summary:</b>	May act as a codon-independent translation release factor that has lost all stop codon specificity and directs the termination of translation in mitochondrion. May help rescuing stalled mitoribosomes during translation (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.