

Product datasheet for **MC207749**

Rnf121 (NM_029211) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rnf121 (NM_029211) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rnf121
Synonyms: 4930544L10Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207749 representing NM_029211
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGGCAGTGGTGGAGGTGGAGGTTGGAGGAGGTGCTCTTGCAGAAAGGGAGCTAGATGAGGTTGATA
TGTCAGACCTCTCCCAGAGGAGCAGTGGAGGGTCAACATGCTCGAATGCATGCCAACATCGTGGCCA
CGAAGCCATGCACGCTGAAATGGTCCTCATCCTCATTGCTACCTTGGTGGTGGCCAGCTGCTCCTGGTA
CAGTGGAAACAGAGGCATCCCCGTTCTTACAATATGGTGACCCTCTTTCAGATGTGGGTTGTTCCCTCT
ATTTACAGTGAAGCTGCATTGGTGGAGGTTCTTGGTGATCTGGATCTTCTCTCTGCTGTCCACAGCCTT
TGTCACCTTCGAGCCACAAGGAAGCCTCTAGTACAGACAACCCCAAGGTTGGTTTATAAATGGTTTCTG
CTGATCTATAAAATCAGCTATGCCACTGGCATTGTTGGCTACATGGCTGTGATGTTTACCCTCTTTGGTC
TTAACTTATTATTCAAGATTAACCAGAAGATGCAATGGACTTTGGCATTCTCTCTCTCTATGGCCT
CTACTATGGAGTCTTTGAGCGTGACTTTGCAGAAATGTGTGCAGACTATATGGCATCAACCATAGGGTTC
TACAGTGAGTCGGGCATGCCTACCAAGCACCTCTCAGACAGCGTGTGTGCTGTGTGTGGGCAGCAGATCT
TTGTGGATGTCAATGAAGAGGGCATCATCGAGAACACGTACAGGCTTTCCTGCAATCAGTATTCATGA
ATTCTGTATTCGAGGCTGGTGCATCGTGGAAAGAAGCAGACATGTCCCTACTGCAAGGAGAAGGTGGAC
CTGAAGAGGATGTTCAAGCAATCCCTGGAAAGGCCTCATGTCATGTATGGACAGCTGCTAGACTGGCTTC
GATACTTGGTAGCCTGGCAGCCTGTCATCATTGGCCTGGTACAAGGCATCAGCTACATCCTGGGCCTCGA
GTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_029211



[View online »](#)

Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_029211.2</u> , <u>NP_083487.2</u>
RefSeq Size:	1502 bp
RefSeq ORF:	984 bp
Locus ID:	75212
UniProt ID:	<u>Q8R1Z9</u>
Cytogenetics:	7 E2