

Product datasheet for MC207829

Tspan9 (NM_175414) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tspan9 (NM_175414) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tspan9
Synonyms: 6720474K14Rik; 9430079M16Rik; AU018597
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC207829 representing NM_175414

Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAGGGGCTGTCTCTGCTGCTTGAAGTACACGATGTTCTCTTCAATTTGATATTCTGGCTCTGTG
 GCTGTGGGCTGCTTGGAGTGGGCATCTGGCTCTCTGTGTCCCAAGGCAACTTCGCCACCTTCTCCCCAG
 CTTCCCCTCACTGTCCGCAGCCAACCTCGTCATCGCCATCGGCACCATCGTCATGGTGACGGGCTTCCTG
 GGCTGCCTGGGGCCATCAAGGAGAACAAGTGTCTGCTTCTCAGCTTCTTCATCGTGCTGTTGATCATCC
 TGCTGGCAGAGTTGATTCTGATCATCTCTTCTTCGTCTACATGGACAAGGTGAACGAGAATGCCAAGCA
 GGACCTGAAAGAGGGGCTGCTGCTCTACACACAGAGAACAACGTGGGGCTCAAGAACGCCTGGAACATC
 ATCCAGGCCGAGATGCGGTGTTGTGGGGTCACTGACTACACAGACTGGTATCCGGTGCTCGGGGAGAACA
 CAGTTCCTGACAGATGCTGCATGGAGAATCCAGGGCTGTGGGCGAACAGCACCACCCCTTGTGGAG
 AACGGGCTGCTACGAAAAAGTGAAGCTGTGGTTTGATGACAACAAGCACGTGCTAGGCACCGTGGGCATG
 TGCATCTCATCATGCAGATCCTGGGAATGGCCTTTCCATGACACTGTTCCAGCACATCCACGGGACAG
 GTAAAAAGTACGACGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_175414
Insert Size: 720 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_175414.4, NP_780623.1</u>
RefSeq Size:	3930 bp
RefSeq ORF:	720 bp
Locus ID:	109246
UniProt ID:	<u>Q8BJU2</u>
Cytogenetics:	6 F3