

Product datasheet for **MC207866**

Sprn (NM_183147) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sprn (NM_183147) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sprn
Synonyms: C630041J07
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC207866 representing NM_183147
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGA**A**CTGGACTGCTGCCACGTGCTGGGCTCTGCTGCTGGCCCGCCCTTCTCTGTGACAGCTGTTCCG
CCAAGGGCGGCCGCGGAGGCGCTCGAGGCAGTGCCCGGGGGTGC**G**CGGAGGTGCACGCGGGGCATCAAG
AGTACGGTAAGGCCGCGCCCGCTACGGCTCCTCTGCGCGTGGCGCTGCAGGGGCAGCAGCAGGG
GCTGCAGCGGGTGTGGCTGCGGGCCTTGTACCGGCTCTGGCTGGAGGAGGACCTCAGGGCCTGGAGAGC
TAGGCCTGGAGGACGATGAGAATGGGGCAATGGGAGGCAACGGAACCGACCGAGGAGTCTACAGCTACTG
GGCCTGGACTTCCGGCTCAGGGTCCGTGCACAGCCACGTATTTGTCTGCTTCTGGGCGGCACCCTTGGT
GCCCTAGAACTGCTTCCGGCT**TAG**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII

ACCN: NM_183147

Insert Size: 444 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_183147.2</u> , <u>NP_898970.1</u>
RefSeq Size:	3172 bp
RefSeq ORF:	444 bp
Locus ID:	212518
UniProt ID:	<u>Q8BWU1</u>
Cytogenetics:	7 F4
Gene Summary:	Prion-like protein that has PrP(C)-like neuroprotective activity. May act as a modulator for the biological actions of normal and abnormal PrP.[UniProtKB/Swiss-Prot Function]