

## Product datasheet for **MC207868**

### Slc35e3 (NM\_029875) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Slc35e3 (NM\_029875) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Slc35e3  
**Synonyms:** 9330166G04Rik; AI481200  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207868 representing NM\_029875  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGCATCGCTGGCGGACCGGGTCCGGGGCAACGGGCGCATCGCCGCGGGCTCCTGTTCAACCTGCTGG  
 TGTCATCTGCATCGTGTTCCTCAACAAATGGATCTATGTACACCACGGCTTCCCAACATGAGCCTGAC  
 CCTGGTGCATTCGTGGTCACCTGGCTGGGCTTGTACATCTGCCAGAACTGAACATCTTGCCCCAAA  
 AGTCTGCCGCTCTCCAAGCTCCTCCTCGCCCTCAGCTTCTGTGGCTTTGTGGTCTTACCAACCTTT  
 CTCTGCAGAACACACCATAGGCACCTATCAGCTGGCCAAGGCCATGACCACGCCGGTGATCATAGCCAT  
 CCAAACCTTCTGGTACCAGAAGAGATTCTCCGTCAGAATACAGCTCACGCTGATTCTTACTGTAGGT  
 GTAATCCTAAATCTTATTACGATGTGAAGTTTCATTGCTTGGAAATGGTGTGGCCGCCCTTGGTGTGG  
 TAGTCACGTCCTTTATCAAGTGTGGTCCGAGCCAAGCAGCAGCAACTGCAGGTGAACCTCATGCAGCT  
 ACTTTACTACCAGGCCCATGTCTTCTGCCATGCTGCTGGTGGCTGTGCCCTTCTTCGAGCCAGTGT  
 GCAGAGGGAGGAATTTGGTCCCTGGTCAGTTCTGCTTTGCTTATGGTACTGCTGTCTGGAATAATCG  
 CTTTCATGGTGAACCTTATCGATTTACTGGATCATTGGGAACACGTCACCAGTCACCTACAACATGTTCCG  
 ACACTTCAAGTTCTGCATCACTCTCTGCCGAGGATACATTCTATTTAAGGACCCCTGTCGGTAAACCAG  
 GGACTCGGCATCCTATGCACGCTCTTCGGTATCCTCACCTACACCCACTTAAACTCAGCGAGCAGGAAG  
 GCAGTAAGAGTAAGCTGGTCCAACGGCC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_029875  
**Insert Size:** 942 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_029875.2</a> , <a href="#">NP_084151.2</a>
<b>RefSeq Size:</b>	3508 bp
<b>RefSeq ORF:</b>	942 bp
<b>Locus ID:</b>	215436
<b>UniProt ID:</b>	<a href="#">Q6PGC7</a>
<b>Cytogenetics:</b>	10 D2
<b>Gene Summary:</b>	Putative transporter.[UniProtKB/Swiss-Prot Function]