

Product datasheet for MC207924

Gemin8 (NM_146238) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gemin8 (NM_146238) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gemin8
Synonyms:	B130034M20Rik; B930082C09Rik; BC023488; gemin-8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC207924 representing NM_146238 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCTCCAACCTGGAGGGCTAGTGCATCTTGGTATTCTCACCCAGTATATGCAAGATACTGGCAACATT
ATCATCACGCAATGCTTTGGATGCAGGGCCACCAGAATGCCTACAGAAAGTTCAGGGATTCTATTTTAC
ATCCCCGTGGCTCTTCCCTCATGGAGCTTCCCTGGAATCTCCTGCTTATGAGGCTGGGCATCCTTGG
GACTCTCAAGGCCAGCACATGGCGCAGCAGGAGTCTCCCTACCGTGTTCATCCAAAAGCCCTGGGC
AGCCTCTGCGTAATAGCAGTAGAACCCAGGCATCCACAAGAGGGAATGAAGCACGGTGTGAGGAGGAAGA
GCTGGAGTCAGATTCAGATGATGAAGTAGAGTGCACCTGAGCAATATGGAGATCACCGAGGAGCTCCGG
CAGTACTTCGCGCAGACCGAGAGGCACAGAGAGGAGAGAAGGCGACAGCAGCTGGACGACAGCGCC
TGAATGACTATGTGAATGCGGACCATGGCCTGACTTCAACCACCGTAGGTCACTGGAACCCCCATCTGA
GAAGCCCTGGGAGCGCGCCAAGCAGAGATGAAGCGCTTATATGGCAACAGTGCTCCCAAGATCCTGGCC
ATGGAGACTGCTGTACAGCTGAGCTTTGACAAGCACTGTGACAGAAAGCAACCAAGTACTGGCCTGTCA
TTCCCTGAAGTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_146238
Insert Size:	717 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_146238.4 , NP_666350.2
RefSeq Size:	1862 bp
RefSeq ORF:	717 bp
Locus ID:	237221
UniProt ID:	Q8BHE1
Cytogenetics:	X F5
Gene Summary:	<p>The SMN complex plays a catalyst role in the assembly of small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Thereby, plays an important role in the splicing of cellular pre-mRNAs. Most spliceosomal snRNPs contain a common set of Sm proteins SNRPB, SNRPD1, SNRPD2, SNRPD3, SNRPE, SNRPF and SNRPG that assemble in a heptameric protein ring on the Sm site of the small nuclear RNA to form the core snRNP. In the cytosol, the Sm proteins SNRPD1, SNRPD2, SNRPE, SNRPF and SNRPG are trapped in an inactive 6S pICln-Sm complex by the chaperone CLNS1A that controls the assembly of the core snRNP. Dissociation by the SMN complex of CLNS1A from the trapped Sm proteins and their transfer to an SMN-Sm complex triggers the assembly of core snRNPs and their transport to the nucleus (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein (isoform 1).</p>