

## Product datasheet for **MC207984**

### 2410002F23Rik (NM\_025880) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 2410002F23Rik (NM\_025880) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** 2410002F23Rik  
**Synonyms:** MGC27669; MGC69849  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC207984 representing NM\_025880  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCTGAGAGAGTGCTGGTGCCACCCAGATAGGCCGGGGGACCGCTACTACACATACACGGAGTTGT  
 TGGCTATCTCACGGCGTTTCAAGCAGAACCCCAATGAGCTCATGGTCACCTGGATCCTGCGGGTGTATGA  
 CCAGGGAGGCCCTGCCCTGTCCCTGAATTCTGGGGAGCTGGGGCTGCTGGGTGACCTCACCCATGATGCC  
 ATCTTTAACTACCGCTGCAAGGCCCTGCGGGGGCTGGCTGCCAGACTCTCCTGAGCTGGCTGCTGCAGG  
 CCTGGCGTCAGCGTTGGGAATCCTCCCTGCATTTTGGAGCCACCGAGTTACCCCTCAGGCCCTGGACCAC  
 CATGGAGGAAGGCATCCAGTTGGTGCGTGAGCTGGGCATGATTGAGTGGATCTACCTTGACCCAGAAGGG  
 CCCGTGGACCTGGCCCCAGAGGATGTGGCCTTCACTCAAGGCCTGCAGCGGCCCTGCTCACAGCAGCCC  
 CCTCTGAGCTGCGGCTTCACTGGTCAGCCTGCTGGTACGTGGCATGACAGTGTGGAGGCTGTGATGGA  
 GATCCAGACTATTGCTGACGTGGTCTGCTCTGGCGCCAGAGCCATCCAGGCCGACCAAGCTCATGTTG  
 GGGCCCAACCAACTCGTAAGGACCTCTAGGCTGGCTGCTCAGCCACGGTGTACCCCGGAGCAAGTAG  
 ACAGACAGCCCACCAAGTACTCCTAGAACTCTACATCAAAGAAGCCAAGCGCAGTCGCGGCCACCCCAA  
 CTATGGGCTGAACGAGGAGCAGCCCCACCACCCCTTACTCCGACCAGGCTGTGGGGAAGAGCAACCG  
 GTGCGTCATGACT**AA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_025880  
**Insert Size:** 855 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_025880.4</a></u> , <u><a href="#">NP_080156.2</a></u>
<b>RefSeq Size:</b>	2001 bp
<b>RefSeq ORF:</b>	855 bp
<b>Locus ID:</b>	668661
<b>Cytogenetics:</b>	7 B3