

Product datasheet for MC208026

Ier2 (NM_010499) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ier2 (NM_010499) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ier2
Synonyms:	AI317238; Ch1; Pip92
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208026 representing NM_010499 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAAGTACAGAAAGAAGCGCAGCGCATCATGACTCTGTGGTATGGAAGATGTACCACTCTCGCATGC
AGCGAGGTGGCTTGCAGCTCCACCGAGTCTGCAGCTGTCCCTCGTTATGCGCAGCGCTCGAGAGGTCTA
CCTCTCAGCCAAGGTAGAAGCCACCAGCCGAGTTCCCGCCATCCCGCAGGGCTCTTGACCTCGCCTG
CACCCGCCGCGGAAGCCGAAGTTCAGTGGAAGTAGCGTCCCCGAAGCCGTGCAGCTCCGGAGCCCA
TGGATACGCAAGAGGAAGTGCTGCGAGTCCAGGAGACCCCTGCGCTCTGTGACCCGCCCCCGCTAGAGT
CAGCCGCAAGCGCCGAGCAGCAGCGATTTGAGCGACGGTAGTGATGCCGACTGGTACCAAGCAAGAAG
GCCCCGTCTAGAAGAGGTGGAGGGGAGGCGACGTCGGAGGTTCCCGATCGCTGCAGCTTCTCCGGCAC
AAAGCGAAGGTGCCTTCCTAACCTCGCCGCGTCTCCAAAGGCGCTTCTCCAGTCTCCTGAACTGTGG
ACCCGCCGTGCCCGCCGACGCCCCACGTGCGAGGCCAAGCCAGCCTGCCGCCGCGCCGACAATATG
CTCAACGTGCTGGTGCAGCTGTGGTGGCTTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_010499
Insert Size:	666 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_010499.4</u> , <u>NP_034629.3</u>
RefSeq Size:	1541 bp
RefSeq ORF:	666 bp
Locus ID:	15936
UniProt ID:	<u>P17950</u>
Cytogenetics:	8 41.02 cM
Gene Summary:	DNA-binding protein that seems to act as a transcription factor (By similarity). Involved in the regulation of neuronal differentiation, acts upon JNK-signaling pathway activation and plays a role in neurite outgrowth in hippocampal cells (By similarity). May mediate with FIBP FGF-signaling in the establishment of laterality in the embryo (By similarity). Promotes cell motility, seems to stimulate tumor metastasis (By similarity).[UniProtKB/Swiss-Prot Function]