

## Product datasheet for **MC208043**

### St3gal1 (NM\_009177) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	St3gal1 (NM_009177) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	St3gal1
Synonyms:	5330418N22Rik; AI467004; Siat4; Siat4a; St3gal-1; ST3Gall
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208043 representing NM_009177 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGGAGGAAGACCCTCAAGTACCTCACCTTCTCCTGCTCTTCATCTTCCTCACTTCCTTTGCTCTGA  
ATTACTCCAACACGGGAGTACCCAGTGCCTGGTCCCCAAGCAGATGCTCCTGGAATCTCAGAGA  
CCGCAGGTTTCATCAAGTACAGCCGTCACCTGCAGACTGCATCAGCCAGGACAAGGTATCATATTGG  
TTCGACCAGCGCTTCAACAAGACTATGCAGCCCTGCTGACAGTCCACAACGCTCTGATGGAGGAGGACA  
CATACCGTGGTGGCTGAGGCTCCAGCGGAGAGGAAACCAACAACCTGAGCGACACCGTCAAGGAACT  
GTTTCGCCTGGTGCCTGGCAATGTGGATCCTATGTTGAACAAGAGGCTGGTGGGCTGCCGACGCTGTGCA  
GTTGTAGGAAACTCCGGTAACTGAAGGACTCCTCGTATGGGCCTGAGATCGACAGCCATGACTTTGTGC  
TGAGGATGAACAAGGCACCCACAGTGGGTTTTGAGGCAGACGTTGGGAGCCGGACCCACCATCTCGT  
GTATCCTGAGAGCTTCCGGGAGCTGGGAGAGAATGTCAACATGGTCTGGTCCCCTTCAAGACCACCGAC  
CTGCAGTGGGTGATCAGCGCCACCACAGGCACCATCACTCACACCTATGTCCCGTGCCCCGAAGA  
TCAAAGTAAACAAGAGAAGATCCTGATCTACCACCCAGCCTTCATCAAGTATGTCTTTGACAACCTGGCT  
TCAGGGCCATGGGCGGTATCCATCGACTGGCATCCTCTCCATCATCTTCCATTATCTGTGATGAG  
GTGGACTTATATGTTTTGGGGCAGACAGCAAAGGAAATTGGCACCATTACTGGGAGAACAACCCATCAG  
CGGGTGCCTTCCGAAAGACGGGGTTACGATGGTGACTTCGAGTACAACATCACAACCTACCTGGCAGC  
CATCAATAAAATCCGCATCTTCAAGGGAGATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: [https://cdn.origene.com/chromatograms/ja2545\\_e08.zip](https://cdn.origene.com/chromatograms/ja2545_e08.zip)



[View online >](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_009177
<b>Insert Size:</b>	1014 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">BC084730</a> , <a href="#">AAH84730</a>
<b>RefSeq Size:</b>	5784 bp
<b>RefSeq ORF:</b>	1014 bp
<b>Locus ID:</b>	20442
<b>UniProt ID:</b>	<a href="#">P54751</a>
<b>Cytogenetics:</b>	15 D2
<b>Gene Summary:</b>	Responsible for the synthesis of the sequence NeuAc-alpha-2,3-Gal-beta-1,3-GalNAc- found on sugar chains O-linked to Thr or Ser and also as a terminal sequence on certain gangliosides. SIAT4A and SIAT4B sialylate the same acceptor substrates but exhibit different Km values.[UniProtKB/Swiss-Prot Function]