

Product datasheet for **MC208072**

Ndufaf4 (NM_026742) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ndufaf4 (NM_026742) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ndufaf4
Synonyms: 1110007M04Rik; 3000003G13Rik; AW214064
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208072 representing NM_026742
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGGGCTCGCGTGACCCGCGCCTTAAGGAACCTCAACGTGGAGAAGCGAGCGGAGCGAGAGATCAGCA
AGAGGAAGCCCTCCATGGCACCCAAGCACCTTCTACCCGCGACCTCCTGCAGGAGCACCGAGTCAGTA
TCCAGAAATCGAAGAAGTTGTTTCTAAAAAAGATAACAAGCTGCTGTCATTGCTAAGAGATGTATATGTC
GATTCCAAAGATCCGGTGCCTGCCGTTAAAGTTGAACACGACAAGAACCAAAGGAATTCAGAC
TGCCGATAGGAAATCACTTTGATAAGAATATCACGGACATTCCTCAAGGCAAGATTACTGTTGTAGAGGC
ATTGACCTTCTCAATAACCATAAACTCTCTCCAGAAACATGGACTGCTGAGAAAATCGCTCAAGAGTAC
TATTTAGAACTGAAGGATGTAATTCCTTCTCAAATATTTTGTACTTTGAAGTCAAATCTTACCTC
CAGAAGACAGGAAAGCAATACAATCAAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026742
Insert Size: 522 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026742.4</u> , <u>NP_081018.1</u>
RefSeq Size:	3632 bp
RefSeq ORF:	522 bp
Locus ID:	68493
UniProt ID:	<u>Q9D1H6</u>
Cytogenetics:	4 A3
Gene Summary:	May be involved in cell proliferation and survival of hormone-dependent tumor cells. Involved in the assembly of mitochondrial NADH:ubiquinone oxidoreductase complex (complex I) (By similarity).[UniProtKB/Swiss-Prot Function]