

Product datasheet for MC208076

Lysmd3 (NM_030257) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Lysmd3 (NM_030257) Mouse Untagged Clone

Tag: Tag Free

Symbol: Lysmd3

Synonyms: 1110030H10Rik; BC003322

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC208076 representing NM_030257
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCGGGGAGGAATCAGAACCGGACGGTTTCCCTCCCAGGAATTCAGGCCAGTGGCCATGTACTTGCTT
TTGGAATTGTACAGACAATGACATGTTGGAAGAAGATGCTGAAGTCTACGAATTCGGTCCAGGGGAAA
GGAGAAGGTTTGAAGAAGCGCATCGAGAGATAGACTTGATGACATTGTTATTAACCAAGATATCCAG
GAAGGAGACACTCTAAACGCAGTAGCCCTTCAGTACTGCTGTACGGTAGCAGATATCAAAAGAGTTAACA
ATCTCATCAGTGATCAAGACTTTTTTGCTCTTAGGTCTATCAAAATTCAGTTAAAGGTTTAGTTCTTT
GACTGAAACTCTTCATCCTCTGAAAGGAAGACACATTTTACATCCTCCACCTGTTCCGTATTTTCAAGAG
CAAGACATTGTGCCAGCTGATGGTCTCTTTCTCCAGTGAGTCAGCCGGCAGTTTCTAAAAGAAGTGG
ACCGAGACATAGAACAAATAGTCAAGTGTACAGACACCAAGAAGGAGAACCTGAATGAAGTCGTGCTGC
CTTAACAGCACAAACAGGTCCGATTTGAACCTGATAACAAAGCATTACCGAAAGGATCCGTATTACGGA
GCAGACTGGGGGATTGGTTGGTGACAGCTGTAGTGATAATGTTGATAGTGGGTATAATAACACCAAGTGT
TTTATCTGCTGTATTATGAAATTTAGCTAAAGTGGATGTTAGTCACCATTCAACAGTGGGTTCTTCACA
CTTGATCCAGGACTCACACCTCCAACACAGCAGAGAGATGGAAAATGAAATTTGGTCCAACAAAAGGA
ATACCTGTTGGTCAGCAAGATGACCACAACTATATAGGCAAGACCCTCAGGCACATGATGCTCAACACA
AAACGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_030257



[View online »](#)

Insert Size:	918 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_030257.1</u> , <u>NP_084533.1</u>
RefSeq Size:	3139 bp
RefSeq ORF:	918 bp
Locus ID:	80289
UniProt ID:	<u>Q99LE3</u>
Cytogenetics:	13 C3
Gene Summary:	Essential for Golgi structural integrity.[UniProtKB/Swiss-Prot Function]