

Product datasheet for **MC208092**

Fam169b (NM_001013811) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam169b (NM_001013811) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam169b
Synonyms: EG434197
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208092 representing NM_001013811
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCAGTTTATAAAGAGGCTTACCTGTGGACATCCTGGAGGATGACCGGAGGGCTATCAGGCAGCAG
CCGAGGCTTACTATGAGATGCTCAGGGAAGGAGCCAGACCTCAGCAGAGTTATCTCCCTTCCACTGG
GGAGCAGGTCAGACTTGAGACTTCATCTCTTTGTTTTGCACGATATACCGTGATGAACCTCAACATAAA
ATACTGGGTTTGGTTAATCCCCAGGACACAAGACAGTTGTGGCTGTTTACTTAAAGGAATCCTGGTGGT
CCATCGAAGACATCCTAAGAACTTCTGATCCCACGCGGGAAGGGCTCATGAAGTTTCAGTCTTTGGGGA
AAGGATTGTTCTCTTTGTTCTCAACGTCATTGTTTTGGAAGACTGGAGAGAAGGCTGCACATCGATGAC
ATGTTTTTCTTACCTCACCTGCGAAGGAGCAAGCTAAGATTCTGTGAAAGATGGAGCCGCGGTTGCAT
TTTACTCAGTCAAAATGAAAGGCAGCCTCTGTGGTGATGGTACCGGCACATGCTACCTGCTGCCTGTCTT
GGACACCGTGTGTTGTGAGGAGGAAGAACCCTTGCCAAGGCTGGGACAGCCATGCTGAGGGACTTCTGT
GACACCTTCAAGGAGACGAAGCTCTGGGATCAGTTGTCCGATATCTCCTGCAATGTACAGAGTCCCTCA
GGCAGTTCCTGTTGACCTGTCCAGGCGAGAGAGGGCGCTTGTGGGAGGTGGAGCCCCGGGGCCCTGGGG
CCAGCAGCGTGTTAACATCTGGCTCAAGGTTTATCTTCAGGAGCGGAGACTTCAAGATGGGAGCACCGTG
CACCTAAATGTTCTGAAGAGGACACGGACACTCCCGGGCAGGCTTCTCAGGAGGATGGACCCACACAGT
TTAATCATGGAGAAAGCCACAAGGAGTGGGCTGTTGGAGAACCGGAGAGGACACAGAATGGCCGACGCTG
TGCGCAAGTGTGTGAGGAAGCGGTCAGGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001013811



[View online »](#)

Insert Size: 1014 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013811.2](#), [NP_001013833.1](#)

RefSeq Size: 2831 bp

RefSeq ORF: 1014 bp

Locus ID: 434197

UniProt ID: [Q8CHT6](#)

Cytogenetics: 7