

Product datasheet for MC208122

Ak4 (NM_001177605) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ak4 (NM_001177605) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ak4
Synonyms: Ak-3; Ak-4; Ak3; Ak3l1; AK 4; D4ErtD274e
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208122 representing NM_001177605
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCTTCCAACTCCTGCGCGGGTCATCCTCGGGCCGCCGGCTCTGGCAAGGGCACCGTTTGCAGAA
 GGATCGCCAGAACTTTGGCTCCAGCATCTCTCCAGCGCCACTTGTTCGGGAGAACCTCAAGACCGG
 CACGGAAGTTGGTGACGTGGCAAACAGTACCTAGAAAAAGGTCTTTGGTTCCAGATCAGTGATCACA
 CGCCTAATGATGTCAGAACTGGAGACTCGGAGTCCCAGCACTGGCTGTTAGATGGATCCCGAGGACAT
 TAGTACAGGCAGAAAGCCCTGGACGGAATCTGTGATGTGGACCTAGTAATCAGTTTGAATATTCCTTTTGA
 GACACTTAAAGATCGTCTGAGCCGACGGTGGATTACCCTTCTAGCGGGAGGGTCTATAACCTGGACTTC
 AACCCACCTCAAGTGCAGGGGATTGATGACATCACTGGTGGAGCCACTTGTCCAACAGGAAGATGATAAAC
 CTGAAGCAGTTGCTGCCAGGCTAAGACGGTACAAGGATGCGGCAAACCAGTCATCGAACTGTACAAGAG
 CCGTGGAGTGCTCCACCAGTTTCTGGGACGGAGACTAACAGAATCTGGCCTTATGTTTACACACTTTTC
 TCCAACAAGATCACACCTATTCAATCAAAGAAGCCTAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001177605
Insert Size: 672 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177605.1</u> , <u>NP_001171076.1</u>
RefSeq Size:	4927 bp
RefSeq ORF:	672 bp
Locus ID:	11639
UniProt ID:	<u>Q9WUR9</u>
Cytogenetics:	4 46.84 cM
Gene Summary:	Involved in maintaining the homeostasis of cellular nucleotides by catalyzing the interconversion of nucleoside phosphates. Efficiently phosphorylates AMP and dAMP using ATP as phosphate donor, but phosphorylates only AMP when using GTP as phosphate donor. Also displays broad nucleoside diphosphate kinase activity.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1-4 encode the same protein.