

## Product datasheet for **MC208127**

### **Alx4 (NM\_007442) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Alx4 (NM\_007442) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Alx4  
**Synonyms:** Ist  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208127 representing NM\_007442  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAACGCTGAGACTTGGCTCTCTTACTGCGAGTCCCCGGCCGCTGCCATGGACGCCTACTACAGCCGG  
TGTCGCAGAGCCGGAAGGTTTCGTCGCCTTTAGGGGATTCAGGAGGGCACAAGTTCGGTACAACCTTT  
TCTGTCCGGCCGGGCCAAAGGACAGGGATTCGGGGACCCAAGAGCCGGGCCGCTACGGCCGGGGCAG  
CAGGACCTGGCGGCACCCTGGAGAGCAGCTCCGGGGCGCGGGGATCCTTTAACAAGTTCAGCCGAGC  
CGCCACCCACAGCCGCGCCGCTCCACCCGCGCCGCTGCGCATCTCTACTTGCAGAGGGGCGCCTG  
CAAAACGCCCCGGACGGTAGCCTCAAGCTCCAGGAGGGCAGCGCGGCCACAACGCGGCCTTGCAGGTC  
CCCTGCTACGCCAAAGAGAGCAACCTGGGTGAGCCAGAGTTGCCTCCTGACTCTGAACCCGTGGGCATGG  
ACAACAGCTACCTCAGTGTGAAGGAGACAGGGGCCAAGGGCCCCAGGACCGGGCCAGCGCTGAAATACC  
AAGCCCCCTGGAAAAGACCGATTCTGAGAGTAACAAAGGCAAGAAGCGGGCAACCGAACCACTTCACT  
AGCTACCAGCTGGAAGAGCTGGAGAAAGTCTCCAGAAGACACACTACCCTGATGTGTATGCACGGGAGC  
AGCTGGCTATGAGGACCGACCTCACCGAGCCCGTGTGCAGTTTGGTTCCAGAACCGGAGGGCCAAAGTG  
GCGAAAGAGGGAGCGCTTTGGGCAGATGCAGCAGGTCCGGACCCACTTCTCCACGCCTATGAGTTGCC  
CTCCTCACCCGGGCTGAGAACTACGCCAGATTAGAACCATCCTGGATTGGCAACAATGGGGCTGCCT  
CGCCAGTGCCAGCCTGTGTGGTCCCCTGTGACCCGGTACCAGCCTGTATGTCCCCTCACGCTCACCTCC  
TGGCTCCGGGGCCAGCAGCGTCTCCGACTTCTTGAGTGTCTCTGGAGCTGGCAGTCACGTGGGCCAGACA  
CACATGGGCAGCCTGTTTGGAGCAGCTGGCATCAGTCCGGGCCTCAATGGCTACGAGATGAACGGCGAGC  
CAGACCGCAAGACCTCAAGCATCGCTGCCTTGCATGAAGGCCAAGGAGCACAGTGCAGCCATCTCTTG  
GGCCACA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_007442
<b>Insert Size:</b>	1200 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_007442.3</a> , <a href="#">NP_031468.1</a>
<b>RefSeq Size:</b>	5203 bp
<b>RefSeq ORF:</b>	1200 bp
<b>Locus ID:</b>	11695
<b>UniProt ID:</b>	<a href="#">O35137</a>
<b>Cytogenetics:</b>	2 51.62 cM
<b>Gene Summary:</b>	Transcription factor involved in skull and limb development.[UniProtKB/Swiss-Prot Function]