

Product datasheet for MC208199

Bst1 (NM_009763) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bst1 (NM_009763) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Bst1
Synonyms:	114/A10; A530073F09; BP-3; Bp3; Bsta1; CD157; Ly65
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208199 representing NM_009763 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGTCCAGGGAGGCCTGTTGTCTCTGTGGCTGTGGCTATGGCTGTCCTTGCTGACAGTGCTTCTGG
GGCAAGAGCGCGGTGGACAGGGGAGGGTACCACTCCTCACCTCCAGAGCATCTTCTGGGCCGTGCGC
AGAGTACACCACGCTGCTGAGCCTGGGGAACAAAACTGCACAGCCATCTGGGAGCCTCAAGGGGGTG
CTGGACAAGGACCCTGCTCTGTGCTCCCCTCAGACTATGATCTCTTCAACCTCCAGGCACCCCA
TTCTAGGGACAAGTCACTGTTCTGGGAGAATAACCACCTGCTCGTTATGAGCTATGGGAGAACACCCG
CCGCTCGTGGCCCTGTGCGACGTGCTGTATGGCAAAGTTGGCGATTTCTTGAGCTGGTGTGACAGGAA
AATGCCTCTGGACTTGACTACCAGTCTGCCCCACATCAGAAGACTGTGAAAACAACGCCGTGGACTCCT
ATTGAAAAGTGCATCCATGCAGTATCAAGAGACAGTTCAGGGGTGATCAACGTGATGCTGAATGGCTC
GGAGCCAAAAGGAGCCTATCCCACGAGAGGGTTTTTTCAGATTTTGAATCCCCTACTTACAAAAGGAT
AAAGTACAAGAATAGAGATCTGGTTATGCATGACGTTGGAGGGCCCAATGTGGAATCCTGTGGCGAAG
GCAGTGTGAAAATCCTGGAAGACAGATTGGAGGCCCTGGGTTCCAGCACAGCTGTATCAATGACTACCG
ACCAGTGAAGTTCTTAATGTGTGTGGACCACAGCACTACCCCGACTGTATCATGAACTCAGCCTCAGCA
TCTATGCGGAGGGAGAGCGCATCTTTCATGCAATAGGAGACGCCAGCCTTCTCATTCTCTTAGTAG
CTTTGGCTTCAAGTTCTCAGGCGTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1729_e01.zip

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_009763
Insert Size:	936 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_009763.3 , NP_033893.2
RefSeq Size:	2336 bp
RefSeq ORF:	936 bp
Locus ID:	12182
UniProt ID:	Q64277
Cytogenetics:	5 23.84 cM
Gene Summary:	Synthesizes the second messengers cyclic ADP-ribose and nicotinate-adenine dinucleotide phosphate, the former a second messenger that elicits calcium release from intracellular stores. May be involved in pre-B-cell growth.[UniProtKB/Swiss-Prot Function]