

Product datasheet for **MC208326**

Crem (NM_001110855) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Crem (NM_001110855) Mouse Untagged Clone
Tag: Tag Free
Symbol: Crem
Synonyms: IC; ICER; ICERI
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208326 representing NM_001110855
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCAAAAGCCCAACATGGCTGTAAGTGGAGATGAACTGCTGCCACAGGTGACATGCCAATTACCAGA
TCCGAGCTCCTACTACTGCTTTGCCACAAGGTGTGGTATGGCTGCCTCACCAGGAAGCCTGCACAGTCC
CCAGCAACTAGCAGAAGAAGCAACTCGCAAGCGGGAGCTGAGGCTGATGAAAAACAGGGAAGCTGCCCGG
GAGTGTGCAGGAAGAAGAAGAATATGTCAAATGTCTTGAAAAATCGTGTGGCTGTGCTTAAAAATCAAA
ACAAGACCTCATTGAGGAAGCAAGGCCCTCAAAGACCTTTATTGCCATAAAGCAGAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001110855
Insert Size: 342 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001110855.1](#), [NP_001104325.2](#)

RefSeq Size: 2010 bp

RefSeq ORF: 342 bp

Locus ID: 12916

UniProt ID: [P27699](#)

Cytogenetics: 18 A1

Gene Summary: This gene encodes a basic-leucine zipper domain-containing protein that localizes to gene promoters, where it binds to the cyclic AMP response element (CRE). Different protein isoforms encoded by this gene may function as either activators or repressors of transcription. Activity of this gene is important in multiple developmental processes, including spermatogenesis. Mutation of this gene causes male infertility. Alternative splicing and promoter usage result in multiple transcript variants for this gene. [provided by RefSeq, Oct 2012]

Transcript Variant: This variant (11) contains multiple differences in the UTRs and coding region, including lacking multiple exons, compared to variant 1. It represents use of an alternate promoter and contains an alternate 5' coding region. The encoded isoform (11) is shorter and has distinct N- and C- termini, compared to isoform 1.