

Product datasheet for MC208330

Crem (NM_001110850) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Crem (NM_001110850) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Crem
Synonyms:	IC; ICER; ICERI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208330 representing NM_001110850 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAGCAAATGTGGCAGGAAAAAGTATATGAGGACAAATGTAAGGCAAATGACCATGGAACAGTTGAAT
CACAGCAGGATCGAAGTGTAAACAGTTCTGTGGCAGAGCATAGCTCTGCTCATATGCAGACTGGTCAAAT
TTCTGTTCTACTCTAGCTCAGGTAGCAACAATTGCAGAGACAGATGATTCTGCAGACTCAGAAGTAATT
GATTTCGCATAAACGTAGAGAAATCTTTCACGAAGACCCTCATATAGAAAAATACTGAATGAACCTTCCT
CTGATGTGCCTGGTATTCCTCAAGATTGAAGAAGAAAAATCAGAGGAAGAAGGGACACCACCTAACATTGC
TACCATGGCAGTACCAACTAGCATATATCAGACTAGCACGGGGCAATACACTGCCACAGGTGACATGCCA
ACTTACCAGATCCGAGCTCCTACTACTGCTTTGCCACAAGGTGTGGTGTGGCTGCCTCACCAGGAAGCC
TGACAGTCCCCAGCAACTAGCAGAAGAAGCAACTCGAAGCGGGAGCTGAGGCTGATGAAAAACAGGGA
AGCTGCCCCGGGAGTGTGCGAGGAAGAAGAAAGAAATATGTCAAATGTCTTGAAAAATCGTGTGGCTGTGCTT
GAAAAACAAAACAAGACCCTCATTGAGGAACCTCAAGGCCCTCAAAGACCTTTATTGCCATAAAGCAGAGT
AA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001110850
Insert Size:	702 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001110850.2, NP_001104320.1</u>
RefSeq Size:	2404 bp
RefSeq ORF:	702 bp
Locus ID:	12916
UniProt ID:	<u>P27699</u>
Cytogenetics:	18 A1
Gene Summary:	<p>This gene encodes a basic-leucine zipper domain-containing protein that localizes to gene promoters, where it binds to the cyclic AMP response element (CRE). Different protein isoforms encoded by this gene may function as either activators or repressors of transcription. Activity of this gene is important in multiple developmental processes, including spermatogenesis. Mutation of this gene causes male infertility. Alternative splicing and promoter usage result in multiple transcript variants for this gene. [provided by RefSeq, Oct 2012]</p> <p>Transcript Variant: This variant (4, also known as alpha gamma) lacks three alternate in-frame exons and uses an alternate splice site in the 3' coding region, compared to variant 1. The encoded isoform (4) is shorter and has a distinct C-terminus, compared to isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>