

Product datasheet for MC208342

Csf2 (NM_009969) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Csf2 (NM_009969) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Csf2
Synonyms:	CSF; Csfgm; Gm-CSF; GMCSF; MGI-IGM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208342 representing NM_009969 Red=Cloning site Blue=ORF

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTGGCTGCAGAATTTACTTTTCTGGGCATTGTGGTCTACAGCCTCTCAGCACCCACCCGCTCACCCA
 TCACTGTCACCCGGCCTTGAAGCATGTAGAGGCCATCAAAGAAGCCCTGAACCTCCTGGATGACATGCC
 TGTCACATTGAATGAAGAGGTAGAAGTCGTCTCTAACGAGTTCTCCTTCAAGAAGCTAACATGTGTGCAG
 ACCCGCCTGAAGATATTCGAGCAGGGTCTACGGGGCAATTTACCAAACCTCAAGGGCGCCTTGAACATGA
 CAGCCAGCTACTACCAGACATACTGCCCCCACTCCGAAACGGACTGTGAAACACAAGTTACCACCTA
 TGCGGATTTATAGACAGCCTTAAACCTTTCTGACTGATATCCCCTTTGAATGCAAAAACAGTCCAA
 AAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_009969
Insert Size:	426 bp


[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC116880 , AA116881
RefSeq Size:	618 bp
RefSeq ORF:	426 bp
Locus ID:	12981
UniProt ID:	P01587
Cytogenetics:	11 32.13 cM
Gene Summary:	<p>Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.</p> <p>[UniProtKB/Swiss-Prot Function]</p>