

Product datasheet for **MC208413**

Eda (NM_001177941) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eda (NM_001177941) Mouse Untagged Clone
Tag: Tag Free
Symbol: Eda
Synonyms: Ed1; Eda-A1; Eda-A2; EDA1; HED; Ta; tabby; Tnlg7c; XLHED
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208413 representing NM_001177941
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGGCTACCCAGAGGTAGAGCGCAGGGAACCCCTGCCTGCGGCAGCGCCAAGGGAGCGGGCAGCCAGG
GCTGCGGCTGTCGCGGGCCCTGCTCGGGCGGGCAAGGGAACAGCTGCCGGCTCTTCTGGGTTTCTT
TGGCCTCTCGCTGGCCCTCCACCTGCTGACGCTGTGCTGTACCTAGAGTTGCGGTCCGAATTGCGGCGG
GAACGGGGAACCGAGTCCCGCCTCGGTGGCCCGGGTGTCTCTGGCACCTCTGGCACCTAAGCAGCCCTG
GGAGCCTCGACCCGGTGGTCCCATCACCCGCCACCTGGGGCAGCCGTCCTTTCAACAGCAGCCTTTGGA
ACCGGGAGAAGATCCACTCCCCCTGACTCCCAGGACCGGCACCAGATGGCCCTCCTGAATTTCTTCTTT
CCTGATGAAAAGGCATATTCTGAAGAGGAAAGTAGGCGTGTTTCGCCCAATAAGAGAAGCAAAAGTGGT
AAGGAGCAGATGGAAAGAAGGCAGGGCCACCTGGGCCAACGGCCCCCAGGACCTCCAGGACCTCCGGG
ACCCAGGGACCTCCAGGGATCCAGGAATTCCTGGGATCCAGGAACAACCTGTTATGGGACCACCTGGC
CCACCTGGCCCTCCTGGTCTCAAGGACCCCTGGCTCCAAGGACCTTCTGGTGTCTGATAAAACTG
GAACCTCGGGAAAATCAGCCAGCTGTGGTGCATCTGCAGGGCCAAGGGTCAGCAATTCAAGTCAAAAATGG
TGGAGTGTCTCAATGACTGGTCTCGCATCACTATGAACCCTAAGGTGTTTAAACTACATCCCCGAGCGGG
GAGCTGGAGGTAAGTGGTGGACGGCACCTACTTCATCTATAGTCAGGTCTACTACATCAACTTCACTGACT
TTGCCAGTACGAGGTGGTGGTGGATGAGAAGCCCTTCTGCAGTGCACCCGAGCATTGAGACAGGGAA
GACCAACTACAACACTTGTATACTGCAGGCGTGTGCCTCCTCAAGGCCAGGCAGAAAAATCGCCGTGAAG
ATGGTGCACGCTGACATCTATCAATATGAGCAAGCACACCACCTTCTCGGGGCCATCAGGCTGGGCG
AAGCCCTGCATCTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001177941
Insert Size:	1137 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177941.1, NP_001171412.1</u>
RefSeq Size:	4916 bp
RefSeq ORF:	1137 bp
Locus ID:	13607
Cytogenetics:	X 43.59 cM
Gene Summary:	<p>Cytokine which is involved in epithelial-mesenchymal signaling during morphogenesis of ectodermal organs. Functions as a ligand activating the DEATH-domain containing receptors EDAR and EDA2R. Isoform TAA binds only to the receptor EDAR, while isoform TA-A2 binds exclusively to the receptor EDA2R (By similarity). May also play a role in cell adhesion (PubMed:10534613).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (6) lacks an in-frame exon and two in-frame small segments in the CDS, as compared to variant 1. The resulting isoform (6) is shorter but has the same N- and C-termini, as compared to isoform 1.</p>