

## Product datasheet for **MC208443**

### Stx2 (NM\_007941) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Stx2 (NM\_007941) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Stx2  
**Synonyms:** AW538950; Epim; G1-536-1; repro34; Syn-2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC208443 representing NM\_007941  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCGGGACCGGCTGCCCGACCTCACGGCGTGTAGGACAAACGACGATGGAGACTGCTGTCGTCATTG  
TGGAGAAGGATCATTTTCATGGACGGTTTCTCCATCAGGTAGAGGAGATTCGAAGCAGCATAGCCAGGAT  
TGCTCAGCATGTAGAAGACGTGAAGAAGAACCACAGCATCATCCTGTCTGCTCCAAACCCAGAAGGAAAA  
ATAAAAGAAGAGCTGGAGACCTGAACAAAGAGATCAAGAAAACCTGCTAACAGGATCCGGGGCAAGCTGA  
AGTCTATTGAGCAGAGCTGTGATCAGGACGAGAATGGGAACCGAACTTCAGTGGATCTGCGGATACGAAG  
GACCCAGCACTCGGTGCTGTACGGAAGTTTGTGGACGTCATGACAGAATACAATGAAGCGCAGATCCTG  
TTCCGGGAGCGAAGCAAAGGCCGCATCCAGCGCCAGCTGGAGATCACTGGGAGGACCACCACTGACGACG  
AGCTGGAAGAGATGCTGGAGAGCGGGAAGCCGTCATCTTCATCTCGGATATTATATCAGATTCACAAAT  
CACTAGGCAAGCTCTCAATGAGATCGAGTCCC GCCACAAAGACATCATGAAGCTGGAGACCAGCATCCGA  
GAGCTGCACGAGATGTTTCATGGATATGGCCATGTTTGTGCGAGACTCAGGGTGAATGGTCAACAACATCG  
AGAGAAATGTGGTGAACCTCTGTAGATTACGTGGAACATGCCAAGGAAGAGACGAAGAAAGCCATCAAATA  
CCAGAGCAAGGCCAGGCCGAAAAAGTGGATAATTGCTGCTGTGGCGGTGGCTGTCTTCCGCTCTGGCT  
CTAATCATTGGCTTGTGCGTTGGCAA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_007941  
**Insert Size:** 870 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_007941.3</a></u> , <u><a href="#">NP_031967.2</a></u>
<b>RefSeq Size:</b>	2937 bp
<b>RefSeq ORF:</b>	870 bp
<b>Locus ID:</b>	13852
<b>Cytogenetics:</b>	5 67.99 cM