

Product datasheet for **MC208464**

F8a (NM_007978) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: F8a (NM_007978) Mouse Untagged Clone
Tag: Tag Free
Symbol: F8a
Synonyms: AI852759; DXUcsf1; F8a1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208464 representing NM_007978
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCTGCGGGCTCTGCGTCCTCCTTGGGCGCGGTGCCTGGCCAGGCTCAGAGGCTGGGGACTTCTTGG
CACGCTATCGGCAGGTGTCCAACAAGCTCAAGAAGCGCTTCTTGAGGAAGCCGAACGTGGCCGAGGCCGG
GGAGCAGTTCGCCAGCTAGCCCGGAGCTGCGCGCCAGGAGTGCCTGCCTATGCTGCCTGGTCCAG
CTGGCTGTGGCGCGTGCACGAGGCGCTTCCATGGGCCCGGGGAAGCGCTGGCCCTCACAGAGCGGG
CCCGACTTTTCTGCGGCAGGAGTGCACGCGCGCCAACGCTTGGGCTGTCCCGCCGCTACGGGGAGCC
TCTGCAGGCGGCCAGCGCCCTCGGCGCCGCTGTGCGCTTGCACCTCGAGCTCGGCCAGCCTGCCGCC
GCTGCTGCACTGTGCCTGGAAGTGTGCTGCCCTTCGCGCAGTGGGCCAGCCGGCCGCTGCTGCCGGTC
ACTTTCAGCGGGTGCAGCTGCACCTGCCCTGATGCCACTGGCCGCACTGCAAGCGCTTGGTGTATGC
TGCCCTCTGCCAGCTGTTGGCGCGGACTACACAGGCGCCCTGGCGCTATTTACACGCATGCAACGCCTG
GCACGGGAGCATGGGGGCCACCCGGTACAGCAACTCGAGCTGCTGCCGAGCCGCCTTCTGGGCCCCAGC
CACCCCTGTGCGGACCCAGCCGAGACCTGTCTGGGCTCGACCTTGGCCCTCCCGAGCCCGCCGACCA
CGCCCCAGGCTCTGTTGCGCCTTACCTGGCACACTCGGTGCCTTGTGACGTCCTAGTCAGGTGTGAG
GTGTCCCGTGTACTGCTGCTGCTCCTGCAACCACCACCTGCCAAGCTGCTGCCCGAGCATGCCAGA
CCCTGGAGAAGTACTCCTGGGAGGCTTTCGATGGCCATGGCCAGGATACCAGCGCCAGCTTCTGAGGA
GCTGTTTCTGTTATTGCACTCTCTGGTATGGCTGCCAAGAAAAGGACACTGAAGGCATCAAGAAGCTG
CAGGTGGAGATGTGGCCACTGCTAACCGCTGAGCAGAACCACCTCCTCCACTCGTTCTGCAGGAAACCA
TCTCTCCCTGACAGGGTGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_007978
Insert Size:	1146 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_007978.3, NP_032004.2</u>
RefSeq Size:	2490 bp
RefSeq ORF:	1146 bp
Locus ID:	14070
Cytogenetics:	X 37.29 cM