

## Product datasheet for **MC208474**

### Fcnb (NM\_010190) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fcnb (NM\_010190) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fcnb  
**Synonyms:** Fcn2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208474 representing NM\_010190  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCCTGGGATCTGCTGCACTATTTCGTCCTTGACCCTGACTGTCCATGCGGCTGGCACATGCCAGAAC  
TGAAGGTCCTAGATCTGGAAGGCTACAAGCAGCTCACCATCCTTCAAGGTTGCCCTGGCTGCCTGGAGC  
TGCAGGGCCCAAGGGAGAGGCAGGTGCCAAAGGAGATAGAGGAGAGAGTGGCCTTCTGGAATTCCTGGA  
AAAGAAGGACCAACTGGACCCAAAGGAAACCAAGGAGAAAAAGGAATACGTGGAGAAAAAGGAGACTCTG  
GGCCCTCTCAGTCATGTGCTACAGGACCTCGGACCTGCAAGGAGTTGCTACCCAGGGCCACTTTCTCAC  
TGGCTGGTATACCATCTATCTGCCAGACTGCAGGCCCTGACTGTGCTGTGTGACATGGACACAGATGGT  
GGAGGCTGGACCGTCTTCCAGAGGAGGCTTGACGGCTCTGTGGACTTCTTTCGGGACTGGACCTCATACA  
AGAGGGGCTTTGGCAGCCAACTAGGGGAGTTCTGGCTGGGGAATGATAATATCCACGCTCTAACCACCCA  
GGGAACCAGTGAGCTGCGGGTGGATCTTTCAGACTTCGAAGGCAAGCATGACTTTGCCAAGTACAGCTCC  
TCCAGATCCAGGGAGAGGCCGAGAAATACAAGCTTATCCTGGGAACTTCCTTGGCGGTGGTGTGGTG  
ACTCTCTGACACCCATAACAACAGGTTATTCTCCACCAAGACCAAGACAATGACGGCAGTACTCCAG  
CTGTGCCATGGGTTACCATGGAGCCTGGTGGTACTCCCAGTGCCACACTTCCAACCTGAATGGCCCTCTAC  
CTGAGGGTCCCCATAAGAGCTATGCAAATGGTGTCAACTGGAAGTCATGGAGAGGGTACAACCTACAGCT  
GCAAGGTTCTGAGATGAAGGTGCGGCTCATCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_010190  
**Insert Size:** 945 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_010190.1</a> , <a href="#">NP_034320.1</a>
<b>RefSeq Size:</b>	1062 bp
<b>RefSeq ORF:</b>	945 bp
<b>Locus ID:</b>	14134
<b>UniProt ID:</b>	<a href="#">O70497</a>
<b>Cytogenetics:</b>	2 A3
<b>Gene Summary:</b>	May function in innate immunity through activation of the lectin complement pathway. Calcium-dependent and GlcNAc-binding lectin (By similarity).[UniProtKB/Swiss-Prot Function]