

Product datasheet for **MC208494**

Fhl3 (NM_010213) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fhl3 (NM_010213) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fhl3
Synonyms: SLIM2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208494 representing NM_010213
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGAGGCATTTGACTGTGCAAAATGCAACGAGTCCTTGTACGGCCGCAAATACATCCAGACAGACA
GTGGCCCTACTGCGTTCGTGCTATGACAACACCTTCGCCAACACCTGTGCTGAGTGCCAGCAGCTCAT
TGGGCATGATTCAAGGGAAGTCTTCTATGAGGACCGCCACTTCCAGAGGGCTGCTCCGCTGCTGCCGA
TGCCAGCGCTCCCTTGCCGATGAGCCCTCACCTGTCAGGACAGTGAGCTTCTGTGTAATGAGTGTACT
GCACAGCCTTCTCGTCCCAGTGTCTGCCTGCGGGGAGACTGTCATGCCTGGGTCCCAGGAACTGGAGTA
TGGAGGCCAGACGTGGCATGAACACTGCTTTCTGTGCAGTGGCTGTGAGCAGCCGCTGGGCTCCCCTCC
TTCGTGCCTGACAAGGGTGCCCACTACTGCGTGCCTTGTATGAGAACAATTTGCTCCTCGGTGTGCC
GCTGCAGCAAGACGCTGACCCAGGGTGGAGTGACATATCGTGATCAGCCCTGGCACCAGAGTGCCTGGT
CTGCACTGGGTGCAAGACGCCCCCTTGCAAGGAGCAGTTCACATCTCGGGATGATGATCCTTACTGTGTG
GCCTGCTTTGGAGAACTTTGCACCAAGTGCAAGCAGCTGCAAGCGCCCATCACAGGTGGGAGCGGGC
GCGGCGAGGGCGCAGGACTCGGTGGAGGCAAGTATGTGTCCTTGAAGACCGACATTGGCACCACAGCTG
CTTTTCTGTGCCGCTGCTCCACTCCCTGGTGGCCAAAGCTTGTACCAAGACGAGATCAAGTTCTA
TGCCAGGGCTGCAGCAAGCAGGCCCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_010213
Insert Size: 870 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010213.3 , NP_034343.2
RefSeq Size:	1685 bp
RefSeq ORF:	870 bp
Locus ID:	14201
UniProt ID:	Q9R059
Cytogenetics:	4 D2.2
Gene Summary:	Recruited by SOX15 to FOXK1 promoters where it acts as a transcriptional coactivator of FOXK1.[UniProtKB/Swiss-Prot Function]