

## Product datasheet for **MC208499**

### Flt3l (NM\_013520) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Flt3l (NM_013520) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Flt3l
Synonyms:	Flt3lg; Ly72L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NM_013520.3 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACAGTGTGGCGCCAGCCTGGAGCCAAATTCCTCCCTGTTGCTGCTGTTGCTGCTGAGTCCTT  
GCCTGCGGGGACACCTGACTGTTACTTCAGCCACAGTCCCATCTCCTCCAACCTCAAAGTGAAGTTTAG  
AGAGTTGACTGACCACCTGCTTAAAGATTACCCAGTCACTGTGGCCGTCATCTTCAGGACGAGAAGCAC  
TGCAAGGCCTTGTGGAGCCTTCTCCTAGCCAGCGCTGGATAGAGCAACTGAAGACTGTGGCAGGGTCTA  
AGATGCAAACGCTTCTGGAGGACGTCAACACCGAGATACATTTTGTACCTCATGTACCTCCAGCCCCT  
ACCAGAATGTCTGCGATTGTCAGACCAACATCTCCACCTCCTGAAGGACACCTGCACACAGCTGCTT  
GCTCTGAAGCCCTGTATCGGGAAGGCCTGCCAGAATTTCTCTCGGTGCCTGGAGGTGCAGTGCCAGCCGG  
ACTCCTCCACCCTGCTGCCCCAAGGAGTCCCATAGCCCTAGAAGCCACGGAGCTCCAGAGCCTCGGCC  
CAGGCAGCTGTTGCTCCTGCTGCTGCTGCTGCCTCTCACACTGGTGTGCTGGCAGCCCGCTGGGGC  
CTTCGCTGGCAAAGGGCAAGAAGGAGGGGGGAGCTCCACCCTGGGGTGCCCTCCCCTCCCATCCCTAG

**ACGCGT**ACGCGGGCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms:	<a href="https://cdn.origene.com/chromatograms/ja1323_d03.zip">https://cdn.origene.com/chromatograms/ja1323_d03.zip</a>
Restriction Sites:	SgfI-MluI
ACCN:	NM_013520
Insert Size:	699 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_013520.3](#), [NP\\_038548.3](#)

**RefSeq Size:** 1075 bp

**RefSeq ORF:** 699 bp

**Locus ID:** 14256

**UniProt ID:** [P49772](#)

**Cytogenetics:** 7 29.14 cM

**Gene Summary:** Stimulates the proliferation of early hematopoietic cells by activating FLT3. Synergizes well with a number of other colony stimulating factors and interleukins.[UniProtKB/Swiss-Prot Function]  
 Transcript Variant: This variant (1) represents the longer transcript and is the protein-coding variant.