

## Product datasheet for **MC208549**

### Ghrh (NM\_010285) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ghrh (NM\_010285) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ghrh  
**Synonyms:** Ghr; Ghrf; GRF  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208549 representing NM\_010285  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTGCTCTGGGTGCTCTTTGTGATCCTCATCCTCACCAGTGGCTCCCACTGCTCACTGCCCCCTCAC  
CTCCCTCAGGATGCAGCGACACGTAGATGCCATCTTCACCACCAACTACAGGAACTCCTGAGCCAGT  
GTATGCCCGAAAGTGATCCAGGACATCATGAACAAGCAAGGGGAGAGGATCCAGGAACAAGGGCCAGG  
CTCAGCCGCCAGGAAGACAGCATGTGGACAGAGGACAAGCAGATGACCCTGGAGAGCATCTTGCAGGGAT  
TCCAAGGATGAAGCCTTCAGCGGACGCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_010285

**Insert Size:** 312 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_010285.2](#), [NP\\_034415.1](#)

**RefSeq Size:** 534 bp

**RefSeq ORF:** 312 bp

**Locus ID:** 14601

**UniProt ID:** [P16043](#)

**Cytogenetics:** 2 78.25 cM

**Gene Summary:** This gene encodes a hormone that has stimulatory effects on pituitary growth hormone synthesis and release, and somatotrope expansion. The encoded preproprotein undergoes proteolytic processing to generate the mature peptide that is secreted by hypothalamus. Mice lacking the encoded protein are deficient in the growth hormone, live longer and exhibit growth retardation, enhanced insulin sensitivity and increased xenobiotic metabolism. [provided by RefSeq, Jul 2016]  
Transcript Variant: This variant (1) represents the longer transcript and encodes the functional protein. Both variants 1 and 2 encode the same protein.