

Product datasheet for MC208625

H2-Eb1 (NM_010382) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2-Eb1 (NM_010382) Mouse Untagged Clone
Tag: Tag Free
Symbol: H2-Eb1
Synonyms: Eb; H-2Eb; H2Eb; Ia-4; Ia4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208625 representing NM_010382
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTGTGGCTCCCCAGAGTTCCTGTGTGGCAGCTGTGATCCTGTTGCTGACAGTGTGAGCCCTCCAG
 TGGCTTTGGTCAGAGACTCCAGACCATGGTTTTTGGAACTACTGTAATCTGAGTGTCACTTCTACAACGG
 GACGCAGCGCGTGC GGCTTCTGGAAGATACTTCTACAACCTGGAGGAGAACTGCGCTTCGACAGCGAC
 GTGGGCGAGTTCCGCGCGGTGACCGAGCTGGGGCGGCCAGACGCCGAGAAGCTGGAACAGCCAGCCGGAGT
 TCCTGGAGCAAAAGCGGGCCGAGGTGGACACGGTGTGCAGACACAACATGAGATCTCGGATAAATTCCT
 TGTGCGGCGGAGAGTTGAGCCTACGGTACTGTGTACCCACAAAAGACACAGCCCTGGAACACCACAAC
 CTCCTGGTCTGCTCTGTGAGTACTTCTACCCTGGCAACATTGAAGTCAGATGGTTCCGGAATGGCAAGG
 AGGAGAAAACAGGAATTGTGTCCACGGGCTGGTCCGAAATGGAGACTGGACCTTCAGACACTGGTGAT
 GCTGGAGACGGTTCCTCAGAGTGGAGAGGTTTACACCTGCCAGGTGGAGCATCCCAGCCTGACCGACCT
 GTCACGGTCGAGTGAAAGCACAATCCACATCTGCACAGAACAAGATGTTGAGTGGAGTTGGGGCTTCG
 TGCTGGCCTGCTCTTCTGGAGCGGGCTGTTTCATCTACTTCAGGAACCAGAAAGGACAGTCTGGACT
 TCAGCCAACAGGACTCCTGAGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_010382
Insert Size: 795 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010382.2](#), [NP_034512.2](#)

RefSeq Size: 1669 bp

RefSeq ORF: 795 bp

Locus ID: 14969

Cytogenetics: 17 17.98 cM