

Product datasheet for MC208626

H2-K1 (NM_001001892) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	H2-K1 (NM_001001892) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	H2-K1
Synonyms:	H-2K; H-2K(d); H2-K; K-f
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208626 representing NM_001001892 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTACCGTGCACGCTGCTCCTGCTGTTGGCGGCCGCCCTGGCTCCGACTCAGACCCGCGGGCCAC
ACTCGCTGAGGTATTTTCGTCACCGCGTGTCCCGCCCGCCCTCGGGGAGCCCGGTACATGGAAGTCGG
CTACGTGGACGACACGGAGTTTCGTGCGCTTCGACAGCGACCGGAGAATCCGAGATATGAGCCGCGGGCG
CGGTGGATGGAGCAGGAGGGCCCGAGTATTGGGAGCGGAGACACAGAAAGCCAAGGGCAATGAGCAGA
GTTTCCGAGTGGACCTGAGGACCCTGCTCGGCTACTACAACCAGAGCAAGGGCGGCTCTCACACTATTCA
GGTGATCTCTGGCTGTGAAGTGGGTCCGACGGGCGACTCCTCCGCGGTACCAGCAGTACGCTACGAC
GGCTGCGATTACATCGCCCTGAACGAAGACCTGAAAACGTGGACGGCGCGGACATGGCGCGCTGATCA
CCAAACACAAGTGGGAGCAGGCTGGTGAAGCAGAGAGACTCAGGGCTACCTGGAGGGCAGCTGCGTGGA
GTGGCTCCGAGATACCTGAAGAACGGGAACGCGACGCTGCTGCGCACAGATTCCCAAAGGCCATGTG
ACCCATCACAGCAGACCTGAAGATAAAGTACCCTGAGGTGCTGGGCCCTGGGCTTCTACCTGCTGACA
TCACCCTGACCTGGCAGTTGAATGGGAGGAGCTGATCCAGGACATGGAGCTTGTGGAGACCAGGCTGC
AGGGATGGAACCTTCCAGAAGTGGCATCTGTGGTGGTGCCTCTTGGGAAGGAGCAGTATTACACATGC
CATGTGTACCATCAGGGCTGCCTGAGCCCTCACCTGAGATGGGAGCCTCCTCCACTGCTCTCCA
ACATGGCGACCGTTGCTGTTCTGTTGCTTGGAGCTGCAATAGTCACTGGAGCTGTGGTGGCTTTTGT
GATGAAGATGAGAAGGAGAAACACAGGTGAAAAAGGAGGGGACTATGCTCTGGCTCCAGGCTCCAGACC
TCTGATCTGTCTCTCCAGATTGTAAGTATGATGGTTTCATGACCCTCATTCTTAGCGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2009_e06.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001001892
Insert Size:	1110 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001892.2 , NP_001001892.2
RefSeq Size:	1669 bp
RefSeq ORF:	1110 bp
Locus ID:	14972
UniProt ID:	P01901
Cytogenetics:	17 17.98 cM
Gene Summary:	Involved in the presentation of foreign antigens to the immune system.[UniProtKB/Swiss-Prot Function]