

Product datasheet for **MC208635**

H2-Q2 (NM_010392) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	H2-Q2 (NM_010392) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	H2-Q2
Synonyms:	BE136769; ENSMUSG00000079507; Gm11132; gs14-2; H-2Q2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC208635 representing NM_010392 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCGCTGCGAAGGCTGCTCCTGCTGTTGGTGGCCGCCCTGAACTGACTGAGACACGCGCGGGCTCAC
ACTCGCTGCGGTATTTACCACCGCCGTGCCCGCCTCGGGAGCCCGGTTTCATTATCGTCGG
CTACGTGGACGACACGCAGTTTCGTGCGCTTCGACAGCGACGCGGAGAATCCGAGGATGGAGCCGCGGGCG
CCGTGGATGGAGCAGGAGGGGCCGAGTATTGGGAGCGGAACACAGGTCTCCAAGGAAATGAGCAGA
GTTTCCGAGTGAGCCTGGGACCGCACTGAGCTACTACAACCAGAGCAAGGGAGGCTCTCACACTCTCCA
GTGGTTGGTTGGCTGTGACCTGGGGCCAGACGGAAGCCTACTCCGAGGGTATGAGCAGTCTGCCTACGAT
GGCCGCGATTACCTCGCCCTGAATGAGGATCTGATAACGTGGACAGCGCGGACCTGGCAGCACTGAAGA
CCCGAAGCAAGTTGGAGCAGGCTGGTCTTGACAGAGAAGCGCAGGGCCTACCTGGAGGTGATTGCTTGAC
GTGGCTCCGACAGATACCTGGAGCTCGGGAAGGAGACGCTGCTGCACACAGATCCCCAAAGGCGCATGTG
ACCCATCACCCAGTTCTCAAGGTGATGTACCCTGAGGTGCTGGGCCCTGGGCTTCTACCTGCTGACA
TCACCCGACCTGGCAGTTGAATGGGAGGACCTGACCCAGGACATGGAGCTTGTGGAGCACCAGGCCTTC
AGGGGATGGAACCTTCCAGAAGTGGGCATCTGTGATGGTGCCTTTTGGGAGGAGCCTAGATACACATGC
CATGTGGAACATGAGGGCTGCCTGAGCCCTCACCTGAGATGGGAGCCTCCTCCATCCACTGACTCTT
ACATGGTGATCATTGCTGTTCTGTTGCTTGGAGCTGTGATCATCATTGGAGCTGTGGTGGCTTTTGT
GATGAAAAGAGGGAGAAACACAGGTGAAAAGTAAAGACTACGCTCAAGATCCAGGCAGGGACAGCCCC
CAGAGCTCTGATATCTCTCTCTCTAGATTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_010392



[View online >](#)

Insert Size:	1083 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_010392.2</u> , <u>NP_034522.2</u>
RefSeq Size:	1100 bp
RefSeq ORF:	1083 bp
Locus ID:	15013
Cytogenetics:	17 18.64 cM