

Product datasheet for **MC208640**

H2-T24 (NM_008207) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2-T24 (NM_008207) Mouse Untagged Clone
Tag: Tag Free
Symbol: H2-T24
Synonyms: H-2T24
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208640 representing NM_008207
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTGGGCTCTTATCTTCTGGCTTCTGAGTCACCCCCAGGACGGCGGAGCCCGGTCGCACTCTCTGCATT
ACTGTTACTCAGCTGTGACTGAGCCGGGCCGGGAGTTCCCTTCATTTTTCCGACGTGGCTTCTTAGATAA
CCAGCCCTTCATCCACTACGACAGCAGGAGCATGAAGGCAGAGCCTTGTGCTGACTGGCTGAGGGAAAAT
GCACAGTACTTCACTCATGAGACTGAGTTTTCCACCAATCGGATGAAGATTTCCAGTTGAGCCTGAGAA
ACATTCGGAATACTACAACAGCTCTGGTACCCAGAGTCAGAGAGCAGATGGCTTCCGCCAGCAAGCAGG
CCCTCACACCCTCCAGTTCACATACGGCTGTGAGATGCGGTACAACAGGACCACAGGGCATTGGCAGTAT
GGCTACGATGGCAGTGATTACCTGACCTTGGACTTGGGCTCTATGCAATATATAGCAGCCACCTTCATTG
CTGGCTACACCAAGCGCAAGTGGGAAAACAATGAGTACTGGCTGGAGAAAAGAGAAGACCTACTTGGAGAA
GGAGTGCATCCTATGGCTGCAGAGATACTTAACCATGGGAGGAAAGAACCTTACCCGAAGTACCCTCCA
AAGACAACAGTGACTCATCAGTTCAAACCTAAAGAAAACGTCACCCTGAGGTGCTGGGCCCTGGGCTTCT
ACCCTGCTGACATCACCTGACCTGGCAGTTGAATGGAGAGGAGCTGACCCAGGACACGGAGCTTGTGGA
GACGAGGCCCTCAGGGGATGGGACATTCAGAAGTGGGACGCTGTGGTGGTGCCTTCTGGAGAGGAGCAG
AGATACACATGTCATGTGCAACATGAGGCACTAACGCAGCCCCTTGTCTGAAATGGGAGCCTTCTCAGC
TCACCACGCCACGACAGGAGTATAGCTAGGGGAAGCTGTAGTCTCAAGCAACCTTGTGAGTGTCTCT
TGCGTTTCCTTTATTTGGCATAGTACTTGTGTTTGGACTTACCCGCTACAAAACCTCAGGTCAGAAGAAAG
AACTGGCTGCTCCATCTGTTCTGAAGAAGAGAGGCTC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_008207



Insert Size:	1092 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_008207.3, NP_032233.2</u>
RefSeq Size:	2182 bp
RefSeq ORF:	1092 bp
Locus ID:	15042
Cytogenetics:	17 18.84 cM