

## Product datasheet for **MC208651**

### Hccs (NM\_008222) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hccs (NM\_008222) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Hccs  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208651 representing NM\_008222  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGTGCCTCAGCATCTTCTCCTGCTACTGCAGTAAACGCCTCGAATGCTTCAGATGGTCAGCCAGCTT  
CACCACCTTCAGGATGCCGATGCATAAAGGGCAAAGGAAAGGCTGCCAGTACTGCAGCAACATCTGA  
CCTGACCAGTGAGAGCAAAGCACACACTGTGCCTGCTCACCAAGACCGGGCCTATGATTATGTGGAATGC  
CCTGTTACTGGTCTAGAGCTAAGGATAAAGGAGAGCCTCGATCCTTCAAACCTGATGCCACCACCTAATC  
AGACACCATCCCCAGATCAGCCATTTACTTTGTCTACCTCAAGAGAGGAGTCATCAATCCAAGAGCTGA  
CTCTGAGAAAAAGTGGGTTTACCCTTCTGAACAGATGTTCTGGAATGCCATGTTACGAAAAGGGTGAAG  
TGGAAAGATGATGATATTAGTCAGAAAGATATGTACAATATCATTAGAATTCACAATCAGAATAATGAGC  
AAGCGTGAAGGAGATTTAAAAATGGGAAGCCCTTCATGCTCATGAGTGCCTTGTGGCCCATCGTTGGT  
CCGATTTGGAGGTAAAGCAAGAGAGTATCACCAAGAGCAAGAATTCGTTCCCTGGATGGGGTATGAATTG  
CCTTTTGACCGCCATGACTGGATCATCAATCGTTGTGGCACAGAAGTTAGATATGTGATTGATTACTATG  
ATGGTGGTGGGTGAACAAGGAGTACCAGTTCACTATCCTGGATGTGCGTCTGCCTTTGATTTCATTCTC  
TGCAGTATGGGACAGAATGAAGGTTGCTTGGTGGCGCTGGACTTCG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_008222  
**Insert Size:** 819 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_008222.4</a> , <a href="#">NP_032248.3</a>
<b>RefSeq Size:</b>	2383 bp
<b>RefSeq ORF:</b>	819 bp
<b>Locus ID:</b>	15159
<b>UniProt ID:</b>	<a href="#">P53702</a>
<b>Cytogenetics:</b>	X F5
<b>Gene Summary:</b>	Links covalently the heme group to the apoprotein of cytochrome c.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-3 all encode the same protein. Sequence Note:.