

Product datasheet for MC208656

Gml2 (NM_010416) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gml2 (NM_010416) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gml2
Synonyms: 1700057K19Rik; HemT; Hemt1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC208656 representing NM_010416
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGCTCCCCTTCTTCTGTCAATTCTCATGGGTCTCCCATGGGTGGATACCTCTATCAACAATACCA
 GTGGATTGGTATTCGACAACACTACTGGTGGCTGGACAATGCCATTGAGCCTAGATGGAGTCTGTGCT
 TACATGCCATAAGTGCTACATAAGTAACACCTTCTCATGTCCAAAGCTCTCTGAGTGTCCCTCGAATCTC
 AGACGCTGATGACAGTCTCCTTTCGTGTGAACATTTCGGCTGCTCTATGTTCTCAAGGACTGCACCAAGG
 ACTGCACATTTATTTATAGAGAACACGTGCCTCCAGAACTCCCTAGGGTGCTTAAAGATGTCAAGAACTT
 CTACTTTGTCATGTGCTGTTCAAGTATAACTTGCAATGTGGGAGGCCCACTAATTTAGAGAGGGACTTG
 TTGGATGAAACTTCCATAGAGGAAGAGGTAGTTGCTAGAGCTGAGTGTCTAGGGTGGGTCAACCTCTCC
 TGTGTTTTGCCTTAATCCTCTCCAGCATCATACTGAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_010416
Insert Size: 531 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_010416.2</u> , <u>NP_034546.2</u>
RefSeq Size:	803 bp
RefSeq ORF:	531 bp
Locus ID:	15202
Cytogenetics:	15 D3