

Product datasheet for **MC208658**

Hes3 (NM_008237) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hes3 (NM_008237) Mouse Untagged Clone
Tag: Tag Free
Symbol: Hes3
Synonyms: bHLHb43; Hes-3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208658 representing NM_008237
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGAAGAAGCGCCGTGCACGCATCAACGTGTCACTGGAGCAGCTGAGGTCTCTTCTGGAGAGACACT
 ACTCACATCAGATACGGAAACGAAAGCTGGAGAAGGCCGATATCCTGGAGCTGAGTGTTAAGTACATGAG
 AAGCCTCCAGAACTCACTGCAAGGACTTGGCCAGTACCCAGTGGGGTGGACTACCCGTCGGCTCCAA
 GGCGGCTTGCGTGGCGTGAGCCAGAGGCTTCGGCCCGGAGAGGGAGACAGCGGCCTGCGCTGCCCTTGC
 TGCTCCAGCGCAGGGAAGGCAGCACACCGACAGCGCCAACCCACAGGCGACCTCTGTTCTCAACCTTGT
 CCTCCCGCCATCTGGCCCCCAGTCGCGCTGCAGGTGGCTCCCACTCCCCGAGTCCCCGCTCCCTCTC
 CCTGGAGGTCTCCTTGAGTCTTCCACCGATGTTGTGGCACCATCTGCATCAAAGTCCAGGCCGAGA
 GTACGAGACCCGGGTTTCGCGTGTGGCGACCGTGG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_008237
Insert Size: 528 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_008237.4</u> , <u>NP_032263.2</u>
RefSeq Size:	1829 bp
RefSeq ORF:	528 bp
Locus ID:	15207
UniProt ID:	<u>Q61657</u>
Cytogenetics:	4 83.01 cM
Gene Summary:	Transcriptional repressor of genes that require a bHLH protein for their transcription. [UniProtKB/Swiss-Prot Function]