

## Product datasheet for MC208676

### Hoxa1 (NM\_010449) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hoxa1 (NM\_010449) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Hoxa1  
**Synonyms:** ERA1; Hox-1.6  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208676 representing NM\_010449  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGACAATGCAAGAATGAACCTCTTTCTGGAATACCCCATCCTTGGCAGTGGCGACTCTGGACCTGCT  
 CGGCGCGAGCTTACCCCTCTGACCATGGGATTACAACCTTCCAATCCTGCGCGGTAGTCCAAACAGCTG  
 CGGCGCGGACGACCGCTTCTAGTGGGCAGGGGGTGCAGATCAGCTCGCCCCACCACCACCACCACC  
 CACCACCATCACCACCCAGACGGCTACTTACCAGACTTCTGGAAACCTTGGGATTTCTTATCCCACT  
 CGAGTTGTGGTCCAAGCTATGGCGCGCAGAAGCTTCACTGCGCCTTATGGCCCTATGGATTAATCAGGA  
 AGCAGACGTAAGTGGTGGTACCCCGCTGCGCTCCCGCTGTTTACTCTGGAAACCTCTCGACTCCCATG  
 GTCCAGCATCACCACCACCAGGGTTATGCTGGGGGCACAGTGGGCTCGCCCCAGTACATTACCACCT  
 CATATGGACAAGAGCAGCAGACTCTGGCCCTGGCCAGTATAATAACTCCTTATCCCTCTCCACGCCAG  
 CCACCAAGAAGCCTGTCGTTCCCTGCTTACAGAGCGTCTTCCAGCGCAGACCTTTGACTGGATGAAA  
 GTTAAAAGAAACCCTCCAAAACAGGAAAGTTGGAGAGTACGGCTACGTGGGTCAACCAACGCAGTGC  
 GCACCAATTTACCACCAAGCAGCTCACAGAGCTGGAGAAGGAGTTCCACTTCAACAAGTACCTTACAGC  
 AGCGCGCAGGGTGGAGATTGCCGCTCCCTACAGCTCAATGAGACCCAGGTGAAGATCTGGTCCAGAAT  
 CGCCGCATGAAGCAGAAGAAGCGTGAGAAGGAGGGCTCCTGCCATCTCCCTGCCACTCCTCCTGGCA  
 GCGATGAGAAAACGGAAGAATCATCTGAGAAATCTAGCCCTCGCCAGTCCCTTCTCCGGCATCGTC  
 TACCTCAGACACTCTGACTACCTCCACTGA

**ACGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_010449



[View online »](#)

**Insert Size:** 1011 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_010449.4](#), [NP\\_034579.3](#)

**RefSeq Size:** 2474 bp

**RefSeq ORF:** 1011 bp

**Locus ID:** 15394

**Cytogenetics:** 6 25.4 cM