

## Product datasheet for **MC208742**

### Ifna11 (NM\_008333) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ifna11 (NM\_008333) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ifna11  
**Synonyms:** IFN-[a]11  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC208742 representing NM\_008333  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCTAGGCTCTGTGCTTTCCTCATGATCTTGATAGTGATGAGCTACTGGTCAACCTGTTCTCTAGGAT  
GCGATCTGCCTCACACTTATAACCTCAGGAACAAGAGGGCCTTGAAGTCTGGCACAAATGAGGAGGCT  
CACCCCTTTTCCTGCCTAAAGGACAGGAAGGACTTTGGATTCCCCTTGGAGAAGGTGGATGCCAGCAG  
ATCCAGAAGGCTCAATCCATCCCTGTGCTGCGAGATCTTACCCAGCAGATCTTGAACCTCTTCGCATCAA  
AGGACTCATCTGCTGCTTGAATGCAACCTCCTAGACTCATTCTGCAATGACCTCCACCAGCAGCTCAA  
TGACCTGCAAGGTTGCTGATGCAGCAGGTGGGGTGCAGGAATCTCCCCTGACCCAGGAAGACTCCCTG  
CTGGCTGTGAGGATATACTTCCACAGGATCACTGTGTTCTGAGAGAGAAGAAACACAGCCCTGTGCCT  
GGGAGGTGGTCAGAGCAGAAGTCTGGAGAGCCCTGTCTTCTCTGCCAATGTGCTGGGAAGACTGAGAGA  
AGAGAAGGCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_008333  
**Insert Size:** 573 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_008333.2</a></u> , <u><a href="#">NP_032359.2</a></u>
<b>RefSeq Size:</b>	647 bp
<b>RefSeq ORF:</b>	573 bp
<b>Locus ID:</b>	15964
<b>UniProt ID:</b>	<u><a href="#">Q61716</a></u>
<b>Cytogenetics:</b>	4 42.03 cM
<b>Gene Summary:</b>	Has antiviral and antiproliferative activities (PubMed:15254193). Produced by macrophages and stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase (By similarity). During viral infection, mediates antiviral effect, either directly by inducing interferon-stimulated genes, either indirectly through stimulation of natural killer cells enabling them to control viral replication (PubMed:22912583).[UniProtKB/Swiss-Prot Function]