

Product datasheet for MC208831

Klra12 (NM_010646) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klra12 (NM_010646) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Klra12
Synonyms:	Ly49L; Ly49I1; Ly49I2; Ly49I3; Ly49I4; Ly49p/d
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC208831 representing NM_010646 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGAGTGAGCAGGAGGTCACTTTCTCAGCTGTGAGATTCCATAAGTCTTCAGGGTTCAGAACCGGGTGA
 GGCTTGAGGAGACAGGGAAGCCTCAAAAAGCTGGCCTCAGAGTTCCTGGCAGCTCATTGTGATAGCTCT
 CGGAATCCTCATTCCCTTCGACTGGTAATTGTCTCAGTGCTGTGACAAACATTTTCAGAATAGTCAA
 CAAAATCATGAAGTGCAGGAACTCTAACTGCCACGATAAGTGCAGCACCACCACTCAAAGTGACATCA
 ACTTGAAGGATGAAGTGTGAGTTCTACATCTATAGAGTGTAGGCCAGGCAATGATCTTCTGGAATCCCT
 CCACAAGGAACAGAACAGATGGTACAGTGAAACCAAGACTTTTTTCAGATTCCTCAGACACAGGCAGAG
 GGTTTTGAAAAATATTGGTTCTGTTATGGTATAAAATGTTATTATTCGTCATGGACAGAAAAACGTGGA
 GTGGATGTAAACAGACCTGCCAGATTTCCAGCTTATCCCTTCTGAAGATAGACAATGAGGATGAAGTGA
 GTTCCTTAAGCTCCTGGTTCCTCAGACAGTTGCTGGATTGGATTGTCATATGATAATAAGAAGAAAGAT
 TGGGCATGGATTAACAATGGCCCATCTAACTTGCTTGAACACAATGAAATATAATATAAGAGATGGAG
 GATGTATGTTGTTATCTAAACAAGACTAGACAATGATAACTGTGATAAATCATTCATCTGTATTTGTGG
 GAAGAGATTGGATAAATCCCTCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_010646
Insert Size:	798 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_010646.1, NP_034776.1</u>
RefSeq Size:	870 bp
RefSeq ORF:	798 bp
Locus ID:	16630
Cytogenetics:	6