

Product datasheet for **MC208886**

Lhx8 (NM_010713) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lhx8 (NM_010713) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lhx8
Synonyms: L3; Lhx7
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208886 representing NM_010713
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTATTGGAAGAGCGATCAGATGTTTTGTGTGAAGCTGGAGGGAAGGAAGTCCCGGAGCTGGCGGTTCC
CCCGCAGAAAGTGCCTCCGGGCTCATGTCGGAGGAGTGCAGGCGGCTGCAGCCCTGGCGCCGGGAGGAC
CCGCAAAGGCGCTGGGGAAGAAGGACTGGTGAATCCCGAGGGAGCCGGGACGAGGACTCCTGCTCCTCC
TCGGCCCCGCTTTCCCGCTCGTCTCGCCCCAGTCCATGGCCTCGGGCTCCGTGTGCCCGCCAGGCAAGT
GTGTGTGCAGCAGCTGCGGCCTGGAGATTGTGGACAAATACCTCCTCAAGGTGAATGACTTATGCTGGCA
TGTCCGCTGTCTCCTGCAGTGTCTGCAGAACGTCCTGGGAAGGCACACGAGCTGCTACATTAAGGAT
AAAGACATTTTCTGTAAACTCGATTACTTCAGACGGTATGGGACCCGCTGTTCCCGCTGTGGCAGGCACA
TCCACTCGACTGACTGGGTCCGCAGGGCAAGGGCAACGTGTATCACCTGGCCTGCTTCCGCTGCTTTTC
TTGCAAGAGGCAGCTGTCCACGGGAGAGGAGTTCGCCTGGTGGAGGAGAAGGTCTCTGTAGAGTGCAC
TTTGACTGCATGCTGGACAATCTGAAGAGAGAAGTGGAGAACGGTAATGGGATTAGTGTGGAAGGAGCCC
TTCTCACAGAGCAAGACGTCAATCATCCAAAGCCAGCCAAAAGAGCTCGGACCAGCTTCCAGCCGATCA
GCTCCAGGTTATGCAAGCACAGTTCGCTCAGGACAACAACCCAGATGCACAGACCCTCCAGAAACTGGCA
GAAAGGACAGGCTTAAGCAGACGTGTCATACAGGTGTGGTTTTCAGAACTGCAGGGCCCGCCATAAGAAAC
ACGTCACTCCCAACCATTTCTCCTCAGCCCCGTCACAGCTGTCCCTCCTCCAGGCTGTCCCCACCCAT
GTTGGAAGAAATGGCTTATTCTGCCTACGTACCCAGGATGGGACGATGCTGACTGCGCTGCACTCATA
ATGGATGCTCACCAACAACCTCCTGGACTCCAGCCCTTGTACCCCATCCAATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_010713



Insert Size:	1104 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010713.2 , NP_034843.2
RefSeq Size:	1990 bp
RefSeq ORF:	1104 bp
Locus ID:	16875
UniProt ID:	O35652
Cytogenetics:	3 H4
Gene Summary:	Transcription factor involved in differentiation of certain neurons and mesenchymal cells. [UniProtKB/Swiss-Prot Function]