

Product datasheet for MC208904

Lta (NM_010735) Mouse Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Lta (NM_010735) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Lta |
| Synonyms: | h1b382; LT; LT-alpha; LT-[a]; LTalpha; Ltx; LT[a]; TNF-beta; Tnfb; TNFSF1; Tnfsf1b; Tnlg1e |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >MC208904 representing NM_010735 Red =Cloning site Blue =ORF |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGACACTGCTCGGCCGTCTCCACCTCTTGAGGGTGCTTGGCACCCCTCCTGTCTTCCTCCTGGGGCTGC
TGCTGGCCCTGCCTCTAGGGGCCAGGGACTCTCTGGTGCCGTTCTCCGCTGCCAGGACAGCCCATCC
ACTCCCTCAGAAGCACTTGACCCATGGCATCCTGAAACCTGCTGCTCACCTTGTTGGGTACCCAGCAAG
CAGAACTCACTGCTCTGGAGAGCAAGCACGGATCGTGCCTTTCTCCGACATGGCTTCTCTTGAGCAACA
ACTCCCTCCTGATCCCCACCACTGGCCTCTACTTTGTCTACTCCAGGTGGTTTTCTCTGGAGAAAGCTG
CTCCCCCAGGGCCATTCCCACTCCCATCTACCTGGCACACGAGGTCCAGCTCTTTTCTCCCAATACCCC
TTCCATGTGCCTCTCCTCAGTGCAGCAAGTCTGTGTATCCGGGACTTCAAGGACCGTGGGTGCGCTCAA
TGTACCAGGGGGCTGTGTTCTGCTCAGTAAGGGAGACCACTGTCCACCCACACCGACGGCATCTCCCA
TCTACACTTCAGCCCCAGCAGTGATTCTTTGGAGCCTTGTCACTGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

| | |
|--------------------|-----------|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_010735 |
| Insert Size: | 609 bp |


[View online »](#)

| | |
|-------------------------------|--|
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p> |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | BC099464 , AAH99464 |
| RefSeq Size: | 1423 bp |
| RefSeq ORF: | 609 bp |
| Locus ID: | 16992 |
| UniProt ID: | P09225 |
| Cytogenetics: | 17 18.59 cM |
| Gene Summary: | <p>Cytokine that in its homotrimeric form binds to TNFRSF1A/TNFR1, TNFRSF1B/TNFR and TNFRSF14/HVEM (By similarity). In its heterotrimeric form with LTB binds to TNFRSF3/LTBR. Lymphotoxin is produced by lymphocytes and is cytotoxic for a wide range of tumor cells in vitro and in vivo.[UniProtKB/Swiss-Prot Function]</p> |