

Product datasheet for **MC208935**

Mageb2 (NM_031171) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mageb2 (NM_031171) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mageb2
Synonyms: Mage-b2; Mage-rs2; Smage2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208935 representing NM_031171
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCCTAGGGGTCAAAGAGTAAGACCCGCTCCCGTGCAAACGACAGCAGTCACGCAGGGAGGTTCCAG
 TAGTTCAGCCCCTGCAGAGGAAGCAGGGTCTTCTCCTGTTGACCAGAGTGTGGGTCCAGCTTCCCTGG
 TGTTCTGCTCCTCAGGGTGTGAAAACCCCTGGATCTTTTGGTGCAGGTGTATCCTGCACAGGCTCTGGT
 ATAGGTGGTAGAAATGCTGCTGTCTGCCTGATACAAAAAGTTCAGATGGCACCCAGGCAGGACTTCCA
 TTCAGCACACACTGAAAGATCCTATCATGAGGAAGGCTAGTGTGCTGATAGAATTCCTGCTAGATAAGTT
 TAAGATGAAAGAAGCAGTTACAAGGAGTGAATGCTGGCAGTAGTTAACAAGAAGTATAAGGAGCAATTC
 CCTGAGATCCTCAGGAGAACTTCTGCACGCCTAGAATTAGTCTTTGGTCTTGAGTTGAAGGAAATTGATC
 CCAGCACTCATTCTATTTGCTGGTAGGCAAACCTGGGCTTTCCACTGAGGGAAGTTTGAGTAGTAAGT
 GGGGTTGCCTAGGACAGGTCTCCTAATGTCTGCCTAGGTGTGATCTTCATGAAGGTAACCGTGCCACT
 GAGCAAGAGGTCTGGCAATTTCTGCATGGAGTGGGGTATATGCTGGGAAGAAGCACTTGATCTTTGGCG
 AGCCTGAGGAGTTTATAAGAGATGTAGTGCAGGAAAATTACCTGGAGTACCGCCAGGTACCTGGCAGTGA
 TCCCCAAGCTATGAGTTCCTGTGGGGACCCAGAGCCCATGCTGAAACAACCAAGATGAAAGTCTGGAA
 GTTTTAGCTAAAGTCAATGGCACAGTCCCTAGTGCCTTCCCTAATCTCTACCAGTTGGCTCTTAGAGATC
 AGGCAGGAGGGGTGCCAAGAAGGAGAGTTCAAGGCAAGGGTGTTCATTCCAAGGCCCATCCCCAAAGTC
 CTCTAACATGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_031171



[View online »](#)

Insert Size:	993 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031171.1</u> , <u>NP_112448.1</u>
RefSeq Size:	993 bp
RefSeq ORF:	993 bp
Locus ID:	17146
Cytogenetics:	X 40.35 cM