

Product datasheet for **MC208944**

Mc2r (NM_008560) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mc2r (NM_008560) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mc2r
Synonyms: ACTH-R; ACTHR; MC2-R
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC208944 representing NM_008560
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAAGCATATTATCAATTCGTATGAACACACCAATGACACCGCAAGAAATAACTCCGATTGTCCTGATG
TAGTTTTGCCAGAAGAGATATTTTTACAATCTCTGTCATTGGCATATTGGAGAAGTTGATTGTCCTCCT
GGCTGTGATCAAAAATAAAAATCTCCAGTCCCCATGTATTTTTTCATCTGCAGTTGGCCATTTCTGAC
ATGTTGGGAGTCTGTATAAGATCTTGGAAAACATCCTGATCATGTTCCAGAAACATGGGTTATCTTAAGC
CTCGTGGCAGTTTTGAAAGCACAGCAGATGACATCATTGACTGCATGTTTCATCCTCTCTTTGCTGGGCTC
TATCTTACGCCTGTCTGTATTGACAGTACCGTTACATCACCATCTTCCATGCCCTGCAATACCATAGC
ATTGTGACCATGCGCCGACCATCATCACCTAACAATTATCTGGATGTTCTGCACAGGGAGCGGCATCA
CCATGGTGATCTTCTCCACCACATCCCCACAGTGTCTCACCTTACATCGCTGTTCCCTTTGATGCTGGT
TTTTATCTGTGTCTCTACATCCACATGTTCTTACTTGCCCGCTCCCATGCTAGGAAGATCTCTACCCTT
CCTAGAACCAACATGAAGGGTGCCATGACACTAACCATCCTTCTGGAGTCTTTCATCTTGTGGGCC
CCTTTGTGCTCCATGTTCTTAAATGACCTTCTGCCAAATAACCTTACTGTGTTGTACATGTCTCT
CTTCCAGTCAATGGCATGTTGATCATGTGCAATGCAGTTATTGACCCCTTATATATGCCTTTCGGAGC
CCAGAGCTCAGAGATGCATTCAAAGGATGCTCTTCTGCAACCGTATT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_008560
Insert Size: 891 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_008560.3 , NP_032586.1
RefSeq Size:	1639 bp
RefSeq ORF:	891 bp
Locus ID:	17200
UniProt ID:	Q64326
Cytogenetics:	18 41.89 cM
Gene Summary:	Receptor for corticotropin (ACTH). This receptor is mediated by G proteins (G(s)) which activate adenylate cyclase (cAMP).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks two internal exons in the 5' UTR, compared to variant 4. Variants 1, 2, 3 and 4 encode the same protein.