

Product datasheet for **MC209035**

Ninj1 (NM_013610) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Ninj1 (NM_013610) Mouse Untagged Clone

Tag: Tag Free

Symbol: Ninj1

Synonyms: AU024536

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC209035 representing NM_013610
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCGGGCACTGAGGAGTATGAGCTCAACGGCGACCTGCGCCGGGCTCCCCGGTTCCCCGACG
CCTTGCCACCCCGCTGGGGTTTGCGAACCGGCCATCAATGTAAACCATTATGCCAACAAGAAGAGCGC
TGCGGAGAGCATGCTGGACATCGCGCTGCTCATGGCCAACGCGTCGAGCTGAAGCCGTGGTGGAGCAG
GGCAATGATTTGCCTTCTTCGTGCCCTTGTGGTCCTCATCTCTATCTCCCTCGTGTGCAGATAGGAG
TGGCGTGTGCTCATCTTCTGGTCAAGTATGACCTCAACAACCGGCCAAGCACGCCAAGCTGGACTT
TCTTAACAACCTGGCCACGGGACTGGTTTTTCATCATCGTCGTGGTCAACATCTTCATTACGCCTTCGGG
GTCCAGAAGCCTGTAATGGACGTGGCGCCCGGCAG**TAG**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2478_a04.zip

Restriction Sites: SgfI-RsrII

ACCN: NM_013610

Insert Size: 459 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC118937 , AAI18938
RefSeq Size:	1132 bp
RefSeq ORF:	459 bp
Locus ID:	18081
UniProt ID:	O70131
Cytogenetics:	13 A5
Gene Summary:	Homophilic cell adhesion molecule that promotes axonal growth. May play a role in nerve regeneration and in the formation and function of other tissues. Cell adhesion requires divalent cations (By similarity).[UniProtKB/Swiss-Prot Function]