

## Product datasheet for **MC209065**

### Nxph2 (NM\_008752) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nxph2 (NM\_008752) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nxph2  
**Synonyms:** MGC124297; MGC124298  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC209065 representing NM\_008752  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGAGCTTGGCGCCGCTGCCTCTCCTTGTAGTCCCAGGCTTGCTGCAGCTGCTATTTTGTGACAGTGAAG  
 AAGTGATACACAACACAGAGAGTGTGGACTGGGAAGACAGAAGTTCCTGAAACATTGGTTGGAATCT  
 GTTCCATTCAAGGATTACCAGTCTCTGGTCTGTTTGTAAACAGCCTCCAGATCCAAAGCCATCATAT  
 GCAGACAACACGAAAACTTTTGGGATTGGCTGGCAAACATAACAGAAATTCAGGAGCAACTAGCAAGAA  
 CTAAGACAGCCAAATAGTGAAGACAGGAAAAATTTAAAAAATGTTTGGGTGGGGTATTCCACTCTAA  
 CATTAAACTGTCAAACCTAACTTGCTCATCACAGGAAAAATGTTGATCATGGAACCGGAACATTCAGT  
 GTTTATTTCCGACACAATTCAACAGGTCTGGGCAATGTTTCAGTGAGCTTGGTACCTCCCTCTAAAGTAG  
 TAGAATTTGAAATATCCCCCAGTCTACCTGGAGACTAAGGAATCCAAATCTTTAACTGCCATATTGA  
 GTATGAAAAACAGATCGGGCAAAGAAGACCGCTTTGTGCAATTCGACCCATCCAAATCTGCTACCAG  
 GAGCAGACTCAGAGTCATGTGCCTGGTTATGTTCCAAACATTTAAGGTCATCTGCATTACATTATCT  
 TTTACAGTGTGATTATAAACTGTACAAAAGGTTTGTCTGACTACAATTACCATAGTGAGACACCATA  
 CTTATCTTTTGGAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_008752  
**Insert Size:** 786 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_008752.2</a></u> , <u><a href="#">NP_032778.1</a></u>
<b>RefSeq Size:</b>	2579 bp
<b>RefSeq ORF:</b>	786 bp
<b>Locus ID:</b>	18232
<b>UniProt ID:</b>	<u><a href="#">Q61199</a></u>
<b>Cytogenetics:</b>	2 A3
<b>Gene Summary:</b>	May be signaling molecules that resemble neuropeptides and that act by binding to alpha-neurexins and possibly other receptors.[UniProtKB/Swiss-Prot Function]