

Product datasheet for MC209077

Olfr15 (NM_008762) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr15 (NM_008762) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Olfr15
Synonyms:	MOR256-17; OR3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209077 representing NM_008762 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGTGGACAGCAACAGCTCCTCTGGGAGCTTCATTCTGATGGGTGTCTCTGACCATCCCCATCTGG
 AGATCATCTTTTTTGTGTCATCCTGGCCTCTTACTTGTTGACGCTGGTTGGAACTTGACCATCATCCT
 GCTTTCGCGCCTTGATGCTCGGCTCCACACCCCATGTACTTCTCCTCAGCAACCTCTCCTCTCTAGAC
 CTTGCCTTACTACCAGTTCAGTCCCTCAGATGCTGAAAAATTTATGGGGCCAGACAAGACAATCAGCT
 ATGGTGGGTGTGTAACCTCACTCTATGTTTTCCTTTGGCTGGGGCTACTGAGTGCATACTGCTCGTGGT
 GATGGCATTGATCGGTATGTGGCAGTTTGTGGGCCCTGCACTACATGACCGTCATGAATCCTCGCCTC
 TGCTGGGGGCTGGCTGCTATTAGCTGGTTGGGTGGCTTAGGCAACTCCGTGATTCACTCAACATTCACTC
 TCCAGCTCCCATTTTGGCGACACCGAAAAGTGGACAACCTCCTGTGTGAGGTACCGCCATGATTAAT
 GGCCTGTGGAGACACAAGTCTCAATGAGGCGGTGCTCAATGGTGTGTTGTACCTTCTTCACTGTGGTCCCA
 GTAAGCGTCATCCTGGTCTCTTACTGCTTCATTGCTCAGGCAGTGATGAAGATCCGCTCTGTGGAGGGAC
 GTCGAAAGGCTTCAATACGTGTGCTCCCACTTGGTGGTAGTGTCTCTTCTATGGCTCTGCGATCTA
 TGGGTATCTGCTCCAGCTAAGAGCAGTAATCAAAGCAAGGAAAATTCATTTCTCTCTACTCTGTG
 GTCACACCCATGGTGAATCCGCTCATCTATACTCTAAGAAACAAAGAAGTGAAGGGGGCCCTGGGAAGAT
 TGCTGGGGAAGGAAGAGGAGCCAGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_008762


[View online »](#)

Insert Size:	939 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_008762.2 , NP_032788.2
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	18312
UniProt ID:	P23275
Cytogenetics:	16 A1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]