

Product datasheet for MC209088

Olfr30 (NM_146878) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr30 (NM_146878) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Olfr30
Synonyms:	MOR2811; MTPCR07
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC209088 representing NM_146878 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGGCAGCGAGCACTGGAAGTACTCCACCGCTGGCTTTGTGCTGACAAGCCTTTTAAACACAGCCAGA
 CACACCTGTTCTGTTCAGCATGGTGATGTTGGTCTACATCCTTGCTATGGCAGGCAACACGGCCATGGT
 CTTTCTGATTTGGATGGACACCCGACTCCACACCCCATGTACTCCTCCTCAGCCAGCTGTCCTTCCTG
 GACATCTTCTTTACATCTGTCACTGTCCCAAGATGATAGTAGGCTTTCTTTTGGTTGGACGAGCATCT
 CTTTCGGAGGTTGTGGGACAGATGTTTTCTTCATGTTCTTGGGGCTGCAGAGTGCCTCCTGCTGGC
 CCTCATGGCTTATGACCGTTACGTGGCTATCTGCAACCTCTGCGATACCCAGTGCTCATGAGCCGCCGG
 GTGTGCCTGCTCATGGTCGTGGCTCCTGGTTGGGAGGGTCCCTCAACGCCTCCATTGAGCTTCTCTGA
 CCCTTCAGTTCCCTACTGTGGATCACGGAAGATCTCCCACTTCTTGTGAGGTGCCCTCGCTGCTGAT
 GCTGGCCTGTGCAGACACAGAAGCCTACAAGCAGGTAATTTGTGACGGGTGTGGTGGTCTCCTGGTG
 CCTATTACCTTCATTACAGCCTCCTATGCCCTCATCCTGGCTGCTGTGCTCCGGATGCACTCTGTGGAGG
 GCGTCAGAAAGCCCTGGCCACATGCTCCTCCACCTGACAGTCGTCATCTCTTCTATGGGCCCTTGT
 CTACACCTACATGTTACCTGCTTCCTATCACTCGCCGGGCCAAGACGATGTAGTATCTGTCTTTACAGC
 GTTCTCACGCCCATGCTGAACCCCGTCATCTACAGCCTCAGGAACAAGGAAGTTACAGGAGCCATGAAGA
 AAGCAATGAGGAGGTGTGGCGTCAATAGAAATGCCT**GA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_146878


[View online »](#)

Insert Size:	948 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146878.2</u> , <u>NP_667089.1</u>
RefSeq Size:	1066 bp
RefSeq ORF:	948 bp
Locus ID:	18329
Cytogenetics:	11 B1.3
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]